




Unique Innovations

- Heart-rate independent temporal resolution of 75 ms
- Tube power: 30 MHU, 7.3 MHU/min
- Focal spot size: Small: 0.7 x 0.7 mm / Large: 0.9 x 1.1 mm
- Generator power: 200 kW
- kV settings: 70/80/100/120/140 kV
- 3D voxel size: 0.24 mm / 0.33 mm
- Bore size: 78 cm
- Max. table load: Up to 307 kg / 676 lbs
- Reconstruction performance:
Up to 50 ips (1 oncology staging exam with 1000 images in up to 20 sec.)
- SAFIRE*
- Dual Energy
- Adaptive Dose Shield for any spiral CT examination
- X-CARE
- Pediatric CT protocols
- Adaptive ECG-Pulsing including MinDose
- Selective Photon Shield
- Flash Spiral with 458 mm/s scan speed
- 1 kW Scan room heat dissipation
- Tube Guard
- Siemens Remote Services
- FAST Planning
- FAST Cardio Wizard
- FAST Scan Assistant
- FAST Adjust
- CARE kV
- CARE Child
- CARE Configurator
- CARE Contrast III
- CARE Profile
- CARE Dashboard



In the event that upgrades require FDA approval, Siemens cannot predict whether or when the FDA will issue its approval. Therefore, if regulatory clearance is obtained and is applicable to this package, it will be made available according to the terms of this offer.

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and is subject to change without prior notice. Some/All of the features and products described herein may not be available in the United States.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features which do not always have to be present in individual cases.

Siemens reserves the right to modify the design, packaging, specifications, and options described herein without prior notice. Please contact your local Siemens sales representative for the most current information.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

Global Business Unit
Siemens AG
Medical Solutions
Computed Tomography &
Radiation Oncology
Siemensstr. 1
DE-91301 Forchheim
Germany
Phone: +49 9191 18 0
Fax: +49 9191 18 9998

Global Siemens Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
www.siemens.com/healthcare

Legal Manufacturer
Siemens AG
Wittelsbacherplatz 2
DE-80333 Muenchen
Germany

Order No. A91CT-06014-49C1-7600 | Printed in Germany | CC CT WS 06112. | © 06.2011, Siemens AG

www.siemens.com/healthcare

www.siemens.com/somatom-definition-flash

CT Cardiology
Benefits –
SOMATOM
Definition
Flash

Multiply Your Potential in CT Cardiology


Generation Flash

The SOMATOM Definition Flash, Siemens' latest high-end scanner, was especially designed to make CT exams much healthier for your patients.

Its core innovation – the revolutionary Flash Spiral – can be summarized in four words: Flash speed. Lowest dose.

Answers for life.

Sub-mSv Heart



Courtesy of Sakakibara Heart Institute / Tokyo, Japan


Discussions regarding dose in cardiac scanning have changed because dose values below 1 mSv are so frequent that sub-mSv heart can be considered routine for Siemens SOMATOM® Definition Flash owners. Even under unfavorable conditions, the patient exposure with the SOMATOM Definition Flash is less than usual for diagnostic cardiac cath.

"[...] prospectively ECG-triggered high-pitch spiral coronary CTA [CT Angiography] provides excellent image quality at a consistent dose below 1.0 mSv."

Achenbach S et al. Coronary computed tomography angiography with a consistent dose below 1 mSv using prospectively electrocardiogram-triggered high-pitch spiral acquisition. Eur Heart J. 2010 Feb;31(3):340-6.

collimation: 128 x 0.6 mm
 spatial resolution: 0.33 mm
 temp resolution: 75 ms
 scan time: 0.29 s
 scan length: 130 mm
 rotation time: 0.28 s
 100/100 kV, 370 mAs/rotation
 eff. dose: 0.97 mSv

All Heart Rates



Courtesy of Sir Run Run Shaw University HongKong / HongKong, China

The SOMATOM Definition Flash's rotation speed of 0.28 s, together with the two X-ray tubes and detectors, allow for an unmatched temporal resolution of down to 75 ms. It reliably freezes cardiac motion so that even patients with high and irregular heart rates above 80 bpm can now get a reliable diagnosis.


"High-pitch protocols lead to [...] optimal contrast enhancement in patients with heart rates below 65 bpm. In patients with higher heart rates [...] or in case of irregular heart rates, prospectively triggered or retrospectively gated protocols will still be the methods of choice."

Sommer WH et al. Feasibility and radiation dose of high-pitch acquisition protocols in patients undergoing dual-source cardiac CT. AJR Am J Roentgenol. 2010 Dec;195(6):1306-12.

collimation: 128 x 0.6 mm
 spatial resolution: 0.33 mm
 temp resolution: 75 msec
 scan time: 4 s
 scan length: 96 mm
 rotation time: 0.28 s
 100/100 kV, 265 mAs/rotation
 CTDIvol: 26.9 mGy
 DLP: 258 mGycm
 eff. dose: 3.6 mSv

no beta blocker administration
 heartrate: 57 - 153 bpm
 average heartrate: 90 bpm

Myocardial Perfusion



Courtesy of Hospital Sírio-Libanês, São Paulo, Brazil


The SOMATOM Definition Flash has everything necessary in cardiac CT imaging and much more. With access to dynamic myocardial stress perfusion imaging, hemodynamic relevance of the stenosis can be evaluated by using syngo VPCT myocardium.

"The ability to obtain accurate perfusion information, in addition to morphologic information from CT coronary angiography imaging, has significant implications [...]."

Ho KT et al. Stress and rest dynamic myocardial perfusion imaging by evaluation of complete time-attenuation curves with dual-source CT. JACC Cardiovasc Imaging. 2010 Aug;3(8):811-20.

collimation: 128 x 0.6 mm
 spatial resolution: 0.33 mm
 scan time: 31 s
 scan length: 66 mm
 rotation time: 0.35 s
 100/100 kV, 176 mAs/rotation
 DLP: 604 mGycm
 CTDIvol: 90,17 mGy
 eff. dose: 8.5mSv

Dose-neutral Dual Energy



Courtesy of MUSC Medical Center / Charleston, USA

With Siemens' unique Dual Source based Dual Energy (DE), it is now possible to accurately quantify iodine content to visualize organ perfusion. The syngo DE Heart PBV software color codes myocardial perfusion, so that both coronary artery morphology and myocardial perfusion can be assessed in a single CT scan.

"Our study results indicate that cardiac DECT [Dual Energy CT] imaging has a high sensitivity [...] and a good specificity [...] for the qualitative assessment of myocardial perfusion."

Meyer M et al. Cost-effectiveness of substituting dual-energy CT for SPECT in the assessment of myocardial perfusion for the workup of coronary artery disease. Eur J Radiol. 2011 Jan 27. [Epub ahead of print]

collimation: 64 x 0.6 mm
 spatial resolution: 0.33 mm
 scan time: 9 s
 scan length: 142 mm
 rotation time: 0.28 s
 100/150 kV, 165/140 mAs/rotation
 heart rate: 103 bpm
 DLP: 339.38 mGycm
 CTDIvol: 23.9 mGy
 eff. dose: 4.75 mSv