

Case Study

Navigation in the Lung

Supported by *syngo* DynaCT Cardiac

www.siemens.com/healthcare

SIEMENS

Navigation in the Lung

Courtesy of Dr. Wolfgang Hohenforst-Schmidt, Institute of Cardiology-Angiography-Pneumology, Klinikum Coburg, Germany

Thoracic Surgery | Pulmonology

Patient History

A 70-year-old-old female.

Diagnosis

The patient was diagnosed with moderately differentiated lung adenomatous carcinoma four years ago and initially treated with chemotherapy. Currently the patient presents with multiple intrapulmonary nodules, metastases in the sternum and multiple osteoplastic metastases in the thoracic spine. Additionally the patient has cerebral metastases.


The patient was admitted for an update to diagnose the response to chemotherapy. In general the patient suffers from moderate dyspnea.

Treatment

A rapid data acquisition was performed with *syngo* DynaCT Cardiac under deep sedation and jet ventilation. In this intra-procedural 3D-dataset the tumor is clearly visible. Navigation through the bronchi is possible in different ways:

- 1) The 3D volume is overlaid with live fluoroscopy. This way, the physician can see live on the movement of the biopsy forceps toward the tumor and make sure the biopsy captures the tissue (figure 1).

Continued on next page ...



1 3D volume of the lung (*syngo* DynaCT Cardiac) with live fluoroscopy overlay (*syngo* iPilot).

Navigation in the Lung

Courtesy of Dr. Wolfgang Hohenforst-Schmidt, Institute of Cardiology-Angiography-Pneumology, Klinikum Coburg, Germany


Thoracic Surgery | Pulmonology

- 2) In addition to option 1, it is possible to manually mark a path through the bronchi as well as the volume of the tumor on the workstation. This path can then again be overlaid and followed when advancing the forceps to the tumor (figure 2).
- 3) Another option is to do a virtual bronchoscopy. The software currently allows navigation up to the 9th debranching of the bronchial tree (figure 3).

Comments

Intra-procedural 3D imaging of the bronchial tree with *syngo* DynaCT Cardiac together with other software applications (*syngo* iPilot, *syngo* iGuide Toolbox and *syngo* FlyThrough) supports biopsies of early-stage lung cancer. The open structure of the C-arm is well suited for interventions and proves to be useful in the field of pneumology.

Continued on next page ...



2 Manually marked path (*syngo* iGuide Toolbox) through the bronchi and marked tumor volume.


Navigation in the Lung

Courtesy of Dr. Wolfgang Hohenforst-Schmidt, Institute of Cardiology-Angiography-Pneumology, Klinikum Coburg, Germany

Thoracic Surgery | Pulmonology

In the above demonstrated case, however, shows an overlay of the tumor in the 3D data set and the real-time fluoroscopy is very precise due to the fact that the whole procedure – 3D data acquisition and bronchoscopy – is done in the same suite at the same time and nearly in the same position of the diaphragm, guaranteed by the apnea under deep sedation (not general anesthesia) and jet ventilation.

We conclude that with *syngo* DynaCT Cardiac rapid onsite navigation in the lung could soon be a routine application in the bronchoscopy suite. *syngo* DynaCT Cardiac has the power to be the central part of a hybrid interventional pneumology/thoracic surgery suite.



3 Virtual bronchoscopy with *syngo* FlyThrough. Clinical images are from different cases.

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice. Some / All of the features and products described herein may not be available in the United States or other countries.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Siemens reserves the right to modify the design, packaging, specifications and options described herein without prior notice. Please contact your local Siemens sales representative for the most current information.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory-new components, we guarantee the quality of these recycled components.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

For product accessories see:
www.siemens.com/medical-accessories

Siemens Healthcare Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

siemens.com/healthcare