

Table of contents

\\Kremser

verschiedenes

div_lib

Siemens_Rektum

localizer_trufi
t2_tse_sag
ep2d_diff_tra
t2_tse_tra
t2_tse_cor
t1_vibe_dixon_cor_CA

\\Kremser\verschiedenes\div_lib\Siemens_Rektum\localizer_trufi

TA: 0:19 PM: ISO Voxel size: 1.3×1.3×5.0 mmPAT: 2 Rel. SNR: 1.00 : tfi

Properties

Prio recon	Off
Load images to viewer	Off
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	7
Dist. factor	20 %
Position	Isocenter
Orientation	Coronal
Phase enc. dir.	R >> L
Slice group	2
Slices	5
Dist. factor	40 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	3
Slices	5
Dist. factor	40 %
Position	L21.0 P0.0 H0.0 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	50 %
FoV read	400 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	3.73 ms
TE	1.87 ms
Averages	1
Filter	Distortion Corr.(2D), Prescan Normalize
Coil elements	BO1-3;SP3-5

Contrast - Common

TR	3.73 ms
TE	1.87 ms
TD	0 ms
Magn. preparation	None
Flip angle	45 deg
Fat suppr.	Fat sat.

Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution - Common

FoV read	400 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
Base resolution	320
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Reference scan mode	Integrated

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	7
Dist. factor	20 %
Position	Isocenter
Orientation	Coronal
Phase enc. dir.	R >> L
Slice group	2
Slices	5
Dist. factor	40 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	3
Slices	5
Dist. factor	40 %
Position	L21.0 P0.0 H0.0 mm
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	400 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	3.73 ms
Multi-slice mode	Sequential
Series	Interleaved

Geometry - AutoAlign

Slice group	1
Slice group	2
Slice group	3
AutoAlign	---
Position	L21.0 P0.0 H0.0 mm

Geometry - AutoAlign

Orientation	Transversal
Phase enc. dir.	A >> P
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Coronal

Geometry - Navigator**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	ISO
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	P >> A
Transversal	H >> F
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.678755 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	3.73 ms
Segments	1

Physio - PACE

Resp. control	Off
---------------	-----

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	On
Dimension	2D
Reordering	Centric
Asymmetric echo	Off
Flow comp.	No
Multi-slice mode	Sequential
Bandwidth	919 Hz/Px

Sequence - Part 2

Segments	1
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slice-sel.

Sequence - Assistant

Mode	Off
------	-----

\\Kremser\verschiedenes\div_lib\Siemens_Rektum\t2_tse_sag

TA: 2:16 PM: ISO Voxel size: 0.7×0.7×4.0 mmPAT: 2 Rel. SNR: 1.00 : tse

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	29
Dist. factor	0 %
Position	L10.7 P0.0 F2.9 mm
Orientation	S > T2.0
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	10 %
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	3960.0 ms
TE	98 ms
Averages	2
Concatenations	2
Filter	Distortion Corr.(2D), Prescan Normalize
Coil elements	BO1-3;SP4-6

Contrast - Common

TR	3960.0 ms
TE	98 ms
TD	0.0 ms
MTC	Off
Magn. preparation	None
Flip angle	160 deg
Fat suppr.	None
Water suppr.	None
Restore magn.	Off

Contrast - Dynamic

Averages	2
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution - Common

FoV read	250 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
Base resolution	384
Phase resolution	83 %
Phase partial Fourier	Off
Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	Self-calibration

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	29
Dist. factor	0 %
Position	L10.7 P0.0 F2.9 mm
Orientation	S > T2.0
Phase enc. dir.	A >> P
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	3960.0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	2

Geometry - AutoAlign

Slice group	1
AutoAlign	---
Position	L10.7 P0.0 F2.9 mm
Orientation	S > T2.0
Phase enc. dir.	A >> P
Initial Position	L10.7 P0.0 F2.9
L	10.7 mm
P	0.0 mm
F	2.9 mm
Initial Rotation	0.00 deg
Initial Orientation	S > T
S > T	2.0
> C	0.0

Geometry - Saturation

Sat. region	1
Thickness	60 mm
Position	L0.0 A97.9 H0.0 mm
Orientation	Coronal
Shape	Standard
Water suppr.	None
Restore magn.	Off
Special sat.	None

Geometry - Navigator

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	F
Table position	3 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	ISO
Table position	F
Table position	3 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	P >> A
Transversal	H >> F
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.678755 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	3960.0 ms
Concatenations	2

Physio - Cardiac

Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	250 mm
FoV phase	100.0 %
Phase resolution	83 %
Trajectory	Cartesian

Physio - PACE

Resp. control	Off
Concatenations	2

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	On
Dimension	2D
Compensate T2 decay	Off
Reduce Motion Sens.	On
Contrasts	1
Flow comp.	No
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	10.9 ms
Bandwidth	200 Hz/Px

Sequence - Part 2

Define	Turbo factor
Echo trains per slice	8
Phase correction	Automatic
Acoustic noise reduction	None
RF pulse type	Low SAR
Gradient mode	Normal
WARP	Off
Red. EC sensitivity	Off
Turbo factor	23

Sequence - Assistant

Mode	Min flip angle
Min flip angle	150 deg
Allowed delay	30 s

\\Kremser\verschiedenes\div_lib\Siemens_Rektum\lep2d_diff_tra

TA: 2:37 PM: ISO Voxel size: 1.0×1.0×4.0 mmPAT: 2 Rel. SNR: 1.00 : epse

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	30
Dist. factor	10 %
Position	L11.6 P7.7 F6.9 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	4700 ms
TE	64.0 ms
Concatenations	1
Filter	Raw filter, Dynamic Field Corr., Distortion Corr.(2D), Prescan Normalize
Coil elements	BO2,3;SP5

Contrast - Common

TR	4700 ms
TE	64.0 ms
MTC	Off
Magn. preparation	None
Fat suppr.	SPAIR
Fat sat. mode	Strong

Contrast - Dynamic

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Delay in TR	0 ms

Resolution - Common

FoV read	250 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
Base resolution	128
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	On

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	38

Resolution - iPAT

Reference scan mode	EPI/separate
---------------------	--------------

Resolution - Filter Image

Distortion Corr.	On
Mode	2D
Prescan Normalize	On
Dynamic Field Corr.	On
Unfiltered images	Off

Resolution - Filter Rawdata

Raw filter	On
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	30
Dist. factor	10 %
Position	L11.6 P7.7 F6.9 mm
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	4700 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
AutoAlign	---
Position	L11.6 P7.7 F6.9 mm
Orientation	Transversal
Phase enc. dir.	A >> P
Initial Position	L11.6 P7.7 F6.9
L	11.6 mm
P	7.7 mm
F	6.9 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Sat. region	1
Thickness	50 mm
Position	L0.0 A103.5 F3.0 mm
Orientation	Coronal
Fat sat. mode	Strong
Special sat.	None

Geometry - Navigator**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	F
Table position	7 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	ISO
------------------	-----

System - Miscellaneous

Table position	F
Table position	7 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	P >> A
Transversal	H >> F
Coil Combine Mode	Adaptive Combine
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	L11.6 P7.7 F6.9 mm
Orientation	Transversal
Rotation	0.00 deg
A >> P	250 mm
R >> L	250 mm
F >> H	132 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.678755 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	3.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	4700 ms
Concatenations	1

Physio - PACE

Resp. control	Off
Concatenations	1

Diff - Neuro

Diffusion mode	4-Scan Trace
Diff. directions	4
Diffusion Scheme	Monopolar
Diff. weightings	3
b-value 1	50 s/mm ²
b-value 2	500 s/mm ²
b-value 3	1000 s/mm ²
b-value 1	1
b-value 2	2
b-value 3	4
Diff. weighted images	Off
Trace weighted images	On
ADC maps	On
FA maps	Off
Mosaic	Off
Tensor	Off
Noise level	10

Diff - Body

Diffusion mode	4-Scan Trace
Diff. directions	4
Diffusion Scheme	Monopolar
Diff. weightings	3
b-value 1	50 s/mm ²
b-value 2	500 s/mm ²
b-value 3	1000 s/mm ²
b-value 1	1
b-value 2	2
b-value 3	4
Diff. weighted images	Off
Trace weighted images	On
ADC maps	On
Exponential ADC Maps	Off
FA maps	Off
Invert Gray Scale	Off
Calculated Image	Off
b-Value >=	0 s/mm ²
Noise level	10

Diff - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D

Sequence - Part 1

Introduction	On
Optimization	Min. TE
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.72 ms
Bandwidth	1562 Hz/Px

Sequence - Part 2

EPI factor	128
RF pulse type	Normal
Gradient mode	Fast

\\Kremser\verschiedenes\div_lib\Siemens_Rektum\t2_tse_tra

TA: 3:44 PM: ISO Voxel size: 0.7×0.7×3.0 mmPAT: 2 Rel. SNR: 1.00 : tse

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	20
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	100 %
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	5180.0 ms
TE	98 ms
Averages	3
Concatenations	1
Filter	Distortion Corr.(2D), Prescan Normalize
Coil elements	BO1-3;SP3-5

Contrast - Common

TR	5180.0 ms
TE	98 ms
MTC	Off
Magn. preparation	None
Flip angle	160 deg
Fat suppr.	None
Water suppr.	None
Restore magn.	Off

Contrast - Dynamic

Averages	3
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution - Common

FoV read	250 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
Base resolution	384
Phase resolution	80 %
Phase partial Fourier	Off
Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	Self-calibration

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	20
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	5180.0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
AutoAlign	---
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	90.00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	None
Water suppr.	None
Restore magn.	Off
Special sat.	None

Geometry - Navigator**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	ISO
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	P >> A
Transversal	H >> F
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.678755 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	5180.0 ms
Concatenations	1

Physio - Cardiac

Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	250 mm
FoV phase	100.0 %
Phase resolution	80 %
Trajectory	Cartesian

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	On
Dimension	2D
Compensate T2 decay	Off
Reduce Motion Sens.	On
Contrasts	1
Flow comp.	No
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	10.9 ms
Bandwidth	200 Hz/Px

Sequence - Part 2

Define	Turbo factor
Echo trains per slice	14
Phase correction	Automatic
Acoustic noise reduction	None
RF pulse type	Low SAR
Gradient mode	Normal
WARP	Off
Red. EC sensitivity	Off
Turbo factor	23

Sequence - Assistant

Mode	Min flip angle
Min flip angle	150 deg
Allowed delay	30 s

\\Kremser\verschiedenes\div_lib\Siemens_Rektum\t2_tse_cor

TA: 3:50 PM: ISO Voxel size: 0.7×0.7×3.5 mmPAT: 2 Rel. SNR: 1.00 : tse

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	19
Dist. factor	10 %
Position	L6.6 P47.3 F5.2 mm
Orientation	C > T34.3 > S1.1
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	70 %
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	3.5 mm
TR	5720.0 ms
TE	93 ms
Averages	3
Concatenations	1
Filter	Distortion Corr.(2D), Prescan Normalize
Coil elements	BO1-3;SP4-6

Contrast - Common

TR	5720.0 ms
TE	93 ms
MTC	Off
Magn. preparation	None
Flip angle	160 deg
Fat suppr.	None
Water suppr.	None
Restore magn.	Off

Contrast - Dynamic

Averages	3
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution - Common

FoV read	250 mm
FoV phase	100.0 %
Slice thickness	3.5 mm
Base resolution	384
Phase resolution	93 %
Phase partial Fourier	Off
Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	Self-calibration

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	19
Dist. factor	10 %
Position	L6.6 P47.3 F5.2 mm
Orientation	C > T34.3 > S1.1
Phase enc. dir.	R >> L
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	3.5 mm
TR	5720.0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
AutoAlign	---
Position	L6.6 P47.3 F5.2 mm
Orientation	C > T34.3 > S1.1
Phase enc. dir.	R >> L
Initial Position	L6.6 P47.3 F5.2
L	6.6 mm
P	47.3 mm
F	5.2 mm
Initial Rotation	0.00 deg
Initial Orientation	C > T
C > T	34.3
> S	1.1

Geometry - Saturation

Fat suppr.	None
Water suppr.	None
Restore magn.	Off
Special sat.	None

Geometry - Navigator**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	F

Geometry - Tim Planning Suite

Table position	5 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	ISO
Table position	F
Table position	5 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	P >> A
Transversal	H >> F
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.678755 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	5720.0 ms
Concatenations	1

Physio - Cardiac

Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	250 mm
FoV phase	100.0 %
Phase resolution	93 %
Trajectory	Cartesian

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Subtract	Off
----------	-----

Inline - Common

Measurements	1
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Sequence - Part 1

Introduction	On
Dimension	2D
Compensate T2 decay	Off
Reduce Motion Sens.	On
Contrasts	1
Flow comp.	No
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	11.7 ms
Bandwidth	200 Hz/Px

Sequence - Part 2

Define	Turbo factor
Echo trains per slice	13
Phase correction	Automatic
Acoustic noise reduction	None
RF pulse type	Low SAR
Gradient mode	Normal
WARP	Off
Red. EC sensitivity	Off
Turbo factor	25

Sequence - Assistant

Mode	Min flip angle
Min flip angle	150 deg
Allowed delay	30 s

\\Kremser\verschiedenes\div_lib\Siemens_Rektum\t1_vibe_dixon_cor_CA

TA: 2:18 PM: ISO Voxel size: 1.0×1.0×1.0 mmPAT: 4 Rel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Dist. factor	20 %
Position	R12.5 P23.2 H47.3 mm
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	80 %
Slice oversampling	23.1 %
Slices per slab	208
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	1.0 mm
TR	7.22 ms
TE 1	2.39 ms
TE 2	4.77 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D), Prescan Normalize
Coil elements	BO1-3;SP5-7

Contrast - Common

TR	7.22 ms
TE 1	2.39 ms
TE 2	4.77 ms
Flip angle	10.0 deg
Fat suppr.	None
Water suppr.	None
Dixon	On
Dixon evaluation	Off

Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution - Common

FoV read	250 mm
FoV phase	100.0 %
Slice thickness	1.0 mm
Base resolution	256
Phase resolution	80 %
Slice resolution	80 %

Resolution - Common

Phase partial Fourier	Off
Slice partial Fourier	Off
Trajectory	Cartesian
View sharing	Off
Interpolation	Off

Resolution - iPAT

PAT mode	CAIPIRINHA
Accel. factor PE	1
Ref. lines PE	24
Accel. factor 3D	4
Ref. lines 3D	24
Reordering Shift 3D	2
Reference scan mode	GRE/separate
CAIPIRINHA mode	Free

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry - Common

Slab group	1
Slabs	1
Dist. factor	20 %
Position	R12.5 P23.2 H47.3 mm
Orientation	Coronal
Phase enc. dir.	R >> L
Slice oversampling	23.1 %
Slices per slab	208
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	1.0 mm
TR	7.22 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	1

Geometry - AutoAlign

Slab group	1
AutoAlign	---
Position	R12.5 P23.2 H47.3 mm
Orientation	Coronal
Phase enc. dir.	R >> L
Initial Position	R12.5 P23.2 H47.3
R	12.5 mm
P	23.2 mm
H	47.3 mm
Initial Rotation	0.00 deg
Initial Orientation	Coronal

Geometry - Saturation

Fat suppr.	None
Water suppr.	None
Dixon	On
Dixon evaluation	Off
Special sat.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	47 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	ISO
Table position	H
Table position	47 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	P >> A
Transversal	H >> F
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	On - AutoCoilSelect

System - Adjustments

B0 Shim mode	Standard Neck
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	R12.5 P23.2 H47.3 mm
Orientation	Coronal
Rotation	0.00 deg
R >> L	250 mm
F >> H	250 mm
A >> P	52 mm
Reset	Off

System - Tx/Rx

Frequency 1H	63.678755 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

View sharing	Off
Flip angle	10.0 deg
Measurements	1
Burn time-to-center	Off
Temporal interpolation	1
3D centric reordering	Off

Inline - Common

Time to center	70.3 s
----------------	--------

Inline - Inline

Subtract	Off
Measurements	1
StdDev	Off
Liver registration	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Soft Tissue

Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
Measurements	1

Inline - Composing

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

Inline - Maplt

Save original images	On
Maplt	None
Flip angle	10.0 deg
Measurements	1
Contrasts	2
TR	7.22 ms
TE 1	2.39 ms
TE 2	4.77 ms

Sequence - Part 1

Introduction	Off
Dimension	3D
Elliptical scanning	Off
Asymmetric echo	Weak
Contrasts	2
Readout mode	Bipolar
Optimization	Opp/In
Multi-slice mode	Sequential
Bandwidth 1	410 Hz/Px
Bandwidth 2	410 Hz/Px

Sequence - Part 2

RF pulse type	Fast
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
Incr. Gradient spoiling	On

Sequence - Assistant

Mode	Off
------	-----