Diagnostic Ultrasound Transducers  
Semi-Critical Transducer Care Addendum  
for Germany

**US**

**Semi-Critical Transducer Care Addendum for Germany (English) ..... 3**

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All Rights Reserved.

May 2018

The following trademarks are owned by Siemens Medical Solutions USA, Inc. (hereinafter "Siemens"):

ACUSON

All other product names are references to third-party products and are trademarks of their respective companies. Siemens includes references to third-party products in the user documentation for informational purposes only. Siemens does not endorse third-party products referenced in the user documentation. Siemens does not assume responsibility for the performance of third-party products.

Siemens reserves the right to change its products and services at any time. In addition, this publication is subject to change without notice.

**Legal Manufacturer**

Siemens Medical Solutions USA, Inc.  
Ultrasound  
685 East Middlefield Road  
Mountain View, CA 94043  
U.S.A.  
Phone: +1-888-826-9702  
siemens-healthineers.com

**CE Declaration**

This product is provided with a CE marking in accordance with the regulations stated in Council Directive 93/42/EEC of June 14, 1993 concerning Medical Devices. The CE marking only applies to medical devices that have been put on the market according to the above referenced Council Directive.

Unauthorized changes to this product are not covered by the CE marking and the related Declaration of Conformity.

**EU Authorized Representative**

Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen  
Germany

**Siemens Healthineers Headquarters**

Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen  
Germany  
Phone: +49 9131 84-0  
siemens-healthineers.com

## **Semi-Critical Transducer Care Addendum for Germany (English)**

This addendum addresses the recent changes to the Kassenärztliche Bundesvereinigung (KBV), Ultrasound Agreement in Germany. In October 2016, KBV revised the Ultrasound Agreement affecting endosonographic transducers (semi-critical category). Approval for reimbursement of ultrasound examinations requires the attending physician to obtain information from the device manufacturer for at least one material-compatible disinfection process confirmed by expert opinion (KBV Ultrasound Agreement, Chapter C, §9 (4)). A transitional period was extended until 31 March 2018.

To fully comply with this revised KBV requirement, the soaking method with Gigasept FF (new), manufactured by Schülke & Mayr GmbH, was validated and confirmed by expert opinion, demonstrating bactericidal, virucidal, and fungicidal efficacy by reproducibly yielding a 4-log reduction rate. This method applies to all semi-critical transducers, endocavity (endovaginal and endorectal), manufactured by Siemens Medical Solutions USA, Inc.

To fully comply with this revised KBV requirement, the wiping method with Tristel Trio Wipes System, manufactured by Tristel Solutions Ltd., was validated and confirmed by expert opinion, demonstrating bactericidal, virucidal, and fungicidal efficacy by reproducibly yielding a 4-log reduction rate. This method applies to all semi-critical transesophageal (TEE) transducers manufactured by Siemens Medical Solutions USA, Inc.

Refer to this addendum and your ultrasound system's user and reference manuals for complete operating instructions, including procedures for cleaning and disinfecting transducers.

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