



Transarterial embolization of ependymoma metastases in the liver

Oncology Interventions

Study Protocol

Courtesy of

PD Dr. Jan Hinrichs, MD,
Institute for Diagnostic and
Interventional Radiology,
Hannover Medical School,
Germany

Supported by

syngo DynaCT

System & Software

ARTIS pheno VE1
with *syngo* Application
Software VD2

Case Description

Patient history

60-year-old female patient with a history of extraneural myxopapillary ependymoma, initially diagnosed in 2009.

Diagnosis

Liver, lung and lymph node metastases since 2015.

Treatment

First session of transarterial embolization of liver metastases with bland embospheres (100–300 µm).

A catheter was placed in the superior mesenteric artery for arteriography and indirect mesentericoportography followed by arteriography of the celiac trunk.

A cone beam CT (CBCT) scan in the arterial phase was carried out with microcatheter placement in the proper hepatic artery. Embolization of the tumor-feeding branch of the right hepatic artery was performed.

After embolization a final DSA showed a decreased perfusion in the embolized liver arteries.

General comments

With the help of the *syngo* DynaCT 3D vessel map, we can ensure the proper position of the catheter for selective embolization. It enhances the identification of aberrant vessels, making it easy for us to verify that there were no vessel branches to the bowel or stomach distal to the embolization position.

Tips and tricks


To ensure optimum *syngo* DynaCT image quality, proper patient breath-hold is vital.

Explain the importance of the patient's cooperation and give them clear and vigorous breath-hold instructions.

Transarterial embolization of ependymoma metastases in the liver

| Acquisition protocol | 5s DCT Body |
|---------------------------|--|
| Injection protocol | |
| Catheter position | Microcatheter in proper hepatic artery |
| Contrast medium (CM) | 300 mg iodine/mL |
| Dilution (CM/Saline): | 60% |
| Injection volume | 20 mL |
| Injection rate | 2.0 mL/s |
| Duration of injection | 10 s |
| X-ray delay | 5 s |
| Power injector used | No |


| Reconstructions | |
|-----------------------|-------------|
| Name | DCT Body |
| VOI size | Full |
| Slice matrix | 512×512 |
| Kernel type | HU |
| Image characteristics | Normal |
| Reconstruction mode | Nat |
| Viewing preset | DynaCT Body |



DSA scene in embolization position



syngo DynaCT VRT visualization



syngo DynaCT imaging MIP 10 mm
Shows tumor and feeding arteries in axial, coronal and sagittal

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone: +49 9131 84-0
siemens-healthineers.com

The statements by Siemens' customers presented here are based on results that were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption), there can be no guarantee that other customers will achieve the same results.

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this case are available throughout the Siemens sales organization worldwide. All rights reserved.