Atellica IM, ADVIA Centaur e Dimension ExL

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Referências:

1. Boeddinghaus J, Twerenbold R, Nestelberger T, et al. Clinical Validation of a Novel High-Sensitivity Cardiac Troponin I Assay for Early Diagnosis of Acute Myocardial Infarction. *Clin Chem* 2018;64:1347-60.
2. Morrow DA, Cannon CP, Jesse RL, et al. National Academy of Clinical Biochemistry Laboratory Medicine Practice Guidelines: clinical characteristics and utilization of biochemical markers in acute coronary syndromes. *Clin Chem* 2007;53:552-74.
3. Myocardial infarction redefined—a consensus document of The Joint European Society of Cardiology/American College of Cardiology Committee for the redefinition of myocardial infarction. *Eur Heart J* 2000;21:1502-13.
4. Apple FS, Collinson PO. Analytical characteristics of high-sensitivity cardiac troponin assays. *Clin Chem* 2012;58:54-61.
5. Hamm CW, Bassand JP, Agewall S, et al. ESC Guidelines for the management of acute coronary syndromes in patients presenting without persistent ST-segment elevation: The Task Force for the management of acute coronary syndromes (ACS) in patients presenting without persistent ST-segment elevation of the European Society of Cardiology (ESC). *Eur Heart J* 2011;32:2999-3054.
6. Thygesen K, Alpert JS, Jaffe AS, et al. Fourth universal definition of myocardial infarction. *Eur Heart J* 2018;00:1-33.
7. Twerenbold R, Boeddighaus J, Mueller C. Update on high-sensitivity cardiac troponin in patients with suspected myocardial infarction. *Eur Heart J Suppl* 2018;(Supplement G):G2-G10.
8. Keller T, Zeller T, Peetz D, et al. Sensitive troponin I assay in early diagnosis of acute myocardial infarction. *N Engl J Med* 2009;361:868-77.
9. Reichlin T, Hochholzer W, Bassetti S, et al. Early diagnosis of myocardial infarction with sensitive cardiac troponin assays. *N Engl J Med* 2009;361:858-67.
10. Wu AH, Lu QA, Todd J, Moecks J, Wians F. Short- and long-term biological variation in cardiac troponin I measured with a high-sensitivity assay: implications for clinical practice. *Clin Chem* 2009;55:52-8.
11. Reichlin T, Irfan A, Twerenbold R, et al. Utility of absolute and relative changes in cardiac troponin concentrations in the early diagnosis of acute myocardial infarction. *Circulation* 2011;124:136-45.
12. Jaffe AS, Moeckel M, Giannitsis E, et al. In search for the Holy Grail: suggestions for studies to define delta changes to diagnose or exclude acute myocardial infarction: a position paper from the study group on biomarkers of the Acute Cardiovascular Care Association. *European heart journal Acute cardiovascular care* 2014;3:313-6.
13. Jaffe AS, Babuin L, Apple FS. Biomarkers in acute cardiac disease: the present and the future. *J Am Coll Cardiol* 2006;48:1-11.
14. Scharnhorst V, Krasznai K, van 't Veer M, Michels RH. Variation of cardiac troponin I and T measured with sensitive assays in emergency department patients with noncardiac chest pain. *Clin Chem* 2012;58:1208-14.
15. Thygesen K, Mair J, Giannitsis E, et al. How to use high-sensitivity cardiac troponins in acute cardiac care. *Eur Heart J* 2012;33:2252-7.
16. Wu AH, Bolger AF, Hollander JE. Growing pains with the use of high-sensitivity cardiac troponin assays. *J Am Coll Cardiol* 2013;62:1250-1.
17. Clinical and Laboratory Standards Institute. Measurement Procedure Comparison and Bias Estimation Using Patient Samples. Approved guideline - Third Edition. CLSI Document EP09-A3, CLSI, Wayne, PA, 2013.
18. Roffi M, Patrono C, Collet JP, et al. 2015 ESC Guidelines for the management of acute coronary syndromes in patients presenting without persistent ST-segment elevation: Task Force for the Management of Acute Coronary Syndromes in Patients Presenting without Persistent ST-Segment Elevation of the European Society of Cardiology (ESC). *Eur Heart J* 2016;37:267-315.
19. Thygesen K, et al. Fourth universal definition of myocardial infarction. *Eur Heart J*. 2018; 00:1-33.
20. Amsterdam EA, et al. 2014 AHA/ACC Guideline for the Management of Patients with non-ST-elevation acute coronary syndromes: a report of the American College of Cardiology/American Heart Association Task Force on Practice Guidelines. *Circulation*. 2014;130:e344-426.

Siemens Healthcare Diagnósticos LTDA

Av. Mutinga, 3800
4° e 5° andares – Pirituba
05110-902 – São Paulo – SP
CAS 0800 55 4838

www.healthcare.siemens.com.br