

SOMATOM go.Top

Clinical Cases

siemens-healthineers.com/somatom-go-top



International version. Not for distribution or use in the U.S.
Courtesy of Radiologie Herrsching & Gilching, Herrsching, Germany

SIEMENS
Healthineers

SOMATOM go.Top

Stand out in advanced CT procedures

In a market characterized by intense competition, more selective patients, and reimbursement cuts, healthcare providers must find ways to leverage technological advancements and secure income and referrals. To keep the business running, it is crucial for CT departments to differentiate themselves and deliver excellent patient-centered care. We want to help you succeed day after day. This is why we developed the SOMATOM® go. platform. As a member of this family, SOMATOM go.Top supports all users to provide the best possible scan for every type of patient – no matter the clinical demands and challenges. The scanner features a unique tablet-based mobile workflow, user guidance with our GO technologies, and exclusive innovations such as Tin Filter low-dose technology.

SOMATOM go.Top is built for personalization of processes and care, allowing every operator to optimally adapt to the individual patient and indication while interacting with patients in a more personalized way than ever before. Produce excellent results for the full clinical spectrum including Dual Energy imaging, and offer what others cannot – for a successful CT business.





Contents



| | |
|------------------------------------|----|
| Cardiovascular Imaging | 7 |
| Whole body and oncological imaging | 21 |
| Head and Neck Imaging | 37 |
| Pediatric Imaging | 47 |
| Acute Care and Trauma | 55 |
| Musculoskeletal Imaging | 67 |
| CT-guided Intervention Guide&GO | 73 |
| How to get there | 77 |

Cardiovascular Imaging

Cardiovascular Imaging

| | |
|---------------------|------------|
| Scan time | 1 s / 6 s |
| Scan length | 147 mm |
| | 110 kV |
| CTDI _{vol} | 10.3 mGy |
| DLP | 140 mGy cm |
| Heart rate | 61 bpm |



- Adaptive Cardio Sequence for dose-effective ECG-synchronized scanning
- Outstanding image quality due to CARE kV and the 10 kV steps
- Automated, zero-click reconstruction of CPRs with Recon&GO



Curved MPRs



Cinematic VRT



0.8 mm MPRs

Cardiovascular Imaging



| | |
|---------------------|-----------|
| Scan time | 6.8 s |
| Scan length | 173 mm |
| | 70 kV |
| CTDI _{vol} | 3.71 mGy |
| DLP | 62 mGy cm |

- Sharp visualization of aneurysm and coronary ostia with 165 ms temporal resolution
- Minimal exposure thanks to low-kV imaging with 70 kV
- Differential diagnosis of ascending aortic aneurysm

Cardiovascular Imaging

| | |
|---------------------|-----------|
| Scan time | 4.5 s |
| Scan length | 123 mm |
| | 80 kV |
| CTDI _{vol} | 6.02 mGy |
| DLP | 74 mGy cm |

- Clear visualization of aortic valve and coronary ostia in systole thanks to high temporal resolution
- Pre-procedural planning for TAVI



0.8 mm MPRs

Cardiovascular Imaging



Angio view



Cinematic VRT



| | |
|---------------------|------------|
| Scan time | 4.7 s |
| Scan length | 623 mm |
| | 80 kV |
| CTDI _{vol} | 1.88 mGy |
| DLP | 118 mGy cm |

- Low dose scan of heart, aorta and vascular system with 80 kV
- Contrast media reduction potential with low-kV imaging which is especially beneficial for elderly patients


Cardiovascular Imaging

| | |
|---------------------|---------------|
| Scan time | 1.4 s / 5.5 s |
| Scan length | 144 mm |
| | 70 kV |
| CTDI _{vol} | 2.77 mGy |
| DLP | 38 mGy cm |
| Heart rate | 65 bpm |

- Adaptive Cardio Sequence in combination with 70 kV for sub-mSv dose
- Inline CPRs of main coronaries to facilitate communication with cardiologist and referring physician



Curved MPRs



1 mm MPR



Cinematic VRT

Cardiovascular Imaging


| | |
|---------------------|------------|
| Scan time | 9 s |
| Scan length | 540 mm |
| | 90 kV |
| CTDI _{vol} | 15.21 mGy |
| DLP | 830 mGy cm |

- Excellent image quality with low-kV imaging even in obese patients
- Occlusion of abdominal aorta

Cardiovascular Imaging

| | |
|---------------------|------------|
| Scan time | 3.7 s |
| Scan length | 488 mm |
| | 100 kV |
| CTDI _{vol} | 4.85 mGy |
| DLP | 240 mGy cm |

- Less dose and high contrast resolution by using unique 10 kV Steps in clinical routine
- Rule-out of endoleak after endovascular aortic repair



3 mm MPRs

Cardiovascular Imaging



MIP



Cinematic VRT



Curved MPR


| | |
|---------------------|------------|
| Scan time | 9.9 s |
| Scan length | 1347 mm |
| | 80 kV |
| CTDI _{vol} | 1.94 mGy |
| DLP | 263 mGy cm |

- Low radiation dose and high contrast resolution with low-kV imaging

Cardiovascular Imaging

| | |
|---------------------|------------|
| Scan time | 24.7 s |
| Scan length | 1146 mm |
| | 90 kV |
| CTDI _{vol} | 7.56 mGy |
| DLP | 880 mGy cm |

- Combine the Stellar detector with 10 kV Steps and CARE kV for detailed visualization of complex vascular diseases



1 mm sagittal and coronal oblique MPRs




Cinematic VRTs

Spectral Imaging with Dual Energy

| | |
|---------------------|------------|
| Scan time | 38 s |
| Scan length | 1354 mm |
| | AuSn120 kV |
| CTDI _{vol} | 6.2 mGy |
| DLP | 818 mGy cm |

- TwinBeam Dual Energy
- Zero-click reconstructions in Recon&GO
- Improved opacification of smaller vessels with Monoenergetic Plus
- 70 keV image equivalent to 120 kV image impression
- Improved visualization with Dual Energy bone removal
- Evaluation of severity of stenosis after bypass



70 keV



40 keV



DE bone removal




Cinematic VRT with bone removal

Whole Body and Oncological Imaging

Whole Body and Oncological Imaging

| | |
|---------------------|----------------------------------|
| Scan time | 8.95 / 4.6 / 4.6 / 8.95 s |
| Scan length | 654 / 320 / 320 / 654 mm |
| | 120 / 90 / 90 / 90 kV |
| CTDI _{vol} | 10.5 / 10.6 / 10.6 / 10.6 mGy |
| DLP | 685 / 341 / 342 / 698 mGy cm |

- Excellent contrast to noise ratio with powerful and fast low-kV imaging
- Contrast media saving potential using 90 kV



1 mm MPR



Non-contrast



Arterial





Late phase




Whole Body and Oncological Imaging

| | |
|---------------------|------------|
| Scan time | 5.9 s |
| Scan length | 368 mm |
| | 90 kV |
| CTDI _{vol} | 9.61 mGy |
| DLP | 365 mGy cm |

- Provide staging even in most challenging patients
- Reduce metal artifacts using iMAR
- Low-kV imaging even in obese patient without compromises on image quality






Without iMAR



With iMAR




Tin Filter topogram



Whole Body and Oncological Imaging

| | |
|---------------------|------------|
| Scan time | 16 s |
| Scan length | 654 mm |
| | 70 kV |
| CTDI _{vol} | 6.37 mGy |
| DLP | 389 mGy cm |

- Outstanding contrast media enhancement due to 70 kV imaging
- Automated results with Recon&GO and inline rib unfolding




1 mm MPRs





1 mm MPRs

R1
R2
R3
R4
R5
R6
R7
R8
R9
R10
R11
R12



L1
L2
L3
L4
L5
L6
L7
L8
L9
L10
L11
L12

Automated results with Recon&GO


Oncological Imaging

| | |
|---------------------|-----------|
| Scan time | 2 s |
| Scan length | 389 mm |
| | Sn110 kV |
| CTDI _{vol} | 0.45 mGy |
| DLP | 15 mGy cm |

- Ultra low-dose lung scan with Tin Filter, available at many kV levels
- Tin Filter Topogram for further dose savings (DLP 0.07 mGy cm)



Tin Filter Topogram




Whole Body and Oncological Imaging

| | |
|---------------------|------------|
| Scan time | 2 s |
| Scan length | 325 mm |
| | 90 kV |
| CTDI _{vol} | 5.6 mGy |
| DLP | 185 mGy cm |



- Excellent image quality and fast scanning even in challenging positioning
- Obese patient with right arm down




Tin Filter topogram



3 mm MPRs



Monoenergetic Plus
40 keV vs. 70 keV at same windowing level



VNC images

Spectral Imaging with Dual Energy

| | |
|---------------------|------------|
| Scan time | 21 s |
| Scan length | 759 mm |
| | AuSn120 kV |
| CTDI _{vol} | 10.2 mGy |
| DLP | 728 mGy cm |

- TwinBeam Dual Energy enables high-contrast and Virtual Non Contrast applications for advanced diagnostic image quality
- Zero-click postprocessing with Recon&GO

Whole Body and Oncological Imaging

| | |
|---------------------|------------------|
| Scan time | 3 / 2 s |
| Scan length | 265 / 232 mm |
| | 110 / 130 |
| CTDI _{vol} | 6.7 / 7 mGy |
| DLP | 181 / 165 mGy cm |

- Advanced diagnostic image quality enabled by high-contrast applications
- Non-contrast and arterial phase



Non-contrast CT




Arterial phase (1 mm)



Cinematic VRT from arterial phase



p.v. contrast media enhanced TBDE CT (mixed)



VNC (TBDE CT)



Monoenergetic 40 keV

Spectral Imaging with Dual Energy

| | |
|---------------------|------------|
| Scan time | 14 s |
| Scan length | 518 mm |
| | AuSn120 kV |
| CTDI _{vol} | 6.1 mGy |
| DLP | 292 mGy cm |

- Dose neutral acquisition with TwinBeam Dual Energy in portal venous phase
- Zero-click postprocessing with inline results in Recon&GO



Spectral Imaging with Dual Energy

| | |
|---------------------|------------|
| Scan time | 11.4 s |
| Scan length | 507 mm |
| | AuSn120 kV |
| CTDI _{vol} | 9.28 mGy |
| DLP | 430 mGy cm |

- TwinBeam Dual Energy
- Visualization and stone characterization in one examination
- Automatic Recon&GO inline results
- Or dedicated evaluation with *syngo*.CT DE Calculi Characterization directly at the scanner



1 mm MPRs



1 mm MPRs


syngo.CT DE Calculi Characterization

Head and Neck Imaging

Head and Neck Imaging


| | |
|---------------------|------------|
| Scan time | 11.9 s |
| Scan length | 176 mm |
| | 120 kV |
| CTDI _{vol} | 43.9 mGy |
| DLP | 830 mGy cm |

- Excellent grey-white matter differentiation thanks to Stellar detector
- Perfectly oriented MPRs directly from the scanner thanks to Recon&GO landmark detection with ALPHA¹



1 mm MPRs

Head and Neck Imaging



| | |
|---------------------|------------|
| Scan time | 10 s |
| Scan length | 287 mm |
| | 80 kV |
| CTDI _{vol} | 10.4 mGy |
| DLP | 297 mGy cm |

- Outstanding image quality even in shoulder regions thanks to High Power 80
- Obese patient

Head and Neck Imaging

| | |
|---------------------|------------|
| Scan time | 5.9 s |
| Scan length | 169 mm |
| | Sn110 kV |
| CTDI _{vol} | 15.1 mGy |
| DLP | 190 mGy cm |

- Great level of detail with 0.6 mm MPRs thanks to Stellar detector
- Optimized dose levels thanks to Tin Filter



Axial and coronal 0.6 mm MPRs

Head and Neck Imaging



Axial 0.6 mm MPRs

| | |
|---------------------|-----------|
| Scan time | 3.94 s |
| Scan length | 93 mm |
| | Sn100 kV |
| CTDI _{vol} | 6.3 mGy |
| DLP | 81 mGy cm |

- Optimal image quality of sinuses in 0.6 mm slices and minimal dose using Tin Filter

Head and Neck Imaging

| | |
|---------------------|------------------|
| Scan time | 8 / 8 s |
| Scan length | 204 / 202 mm |
| | 120 kV |
| CTDI _{vol} | 45 / 42 mGy |
| DLP | 733 / 669 mGy cm |

- Excellent low-contrast performance for differentiation of gray/white matter
- Restaging after lung cancer resection




1 mm MPRs

Spectral Imaging with Dual Energy



70 keV



40 keV



TBDE bone removal


| | |
|---------------------|------------|
| Scan time | 9 s |
| Scan length | 358 mm |
| AuSn120 kV | |
| CTDI _{vol} | 9.29 mGy |
| DLP | 293 mGy cm |

- TwinBeam Dual Energy
- Zero-click reconstructions in Recon&GO
- Improved visualization with Dual Energy bone removal
- Improved opacification of smaller vessels with Monoenergetic Plus


Head and Neck Imaging

| | |
|---------------------|------------------|
| Scan time | 10 / 15 s |
| Scan length | 197 / 43 mm |
| | 120 / 70 kV |
| CTDI _{vol} | 37 / 162 mGy |
| DLP | 738 / 624 mGy cm |


- Boost your stroke assessment by adding functional information to morphology with 4D imaging
- Experience a routine-ready workflow with neuro perfusion tools directly at the scanner
- Penumbra calculation based on mismatch between CBV and CBF, Tmax and rCBF or any other pair of parameters




2 mm MPR



MIP



CBF




Penumbra

Pediatric Imaging

Pediatrics

| | |
|---------------------|-----------|
| Scan time | 2.7 s |
| Scan length | 249 mm |
| | 70 kV |
| CTDI _{vol} | 0.76 mGy |
| DLP | 21 mGy cm |

- 17 week old child
- Fast acquisition for high quality pediatric scanning



3.0 mm MPR



Follow-up after surgery
3.0 mm MPR

Pediatrics


| | |
|---------------------|-----------|
| Scan time | 5.9 s |
| Scan length | 177 mm |
| | 70 kV |
| CTDI _{vol} | 0.95 mGy |
| DLP | 18 mGy cm |

- CARE Child protocols with 70 kV for low dose in pediatrics
- Follow-up scan after splenectomy


Pediatric Imaging

| | |
|---------------------|-----------|
| Scan time | 3.8 s |
| Scan length | 348 mm |
| | 70 kV |
| CTDI _{vol} | 1.58 mGy |
| DLP | 61 mGy cm |

- Experience better iodine enhancement and lower dose thanks to dedicated CARE Child protocols
- 11 year old with laceration of the liver



3.0 mm MPRs



Pediatric Imaging

| | |
|---------------------|------------|
| Scan time | 7.5 s |
| Scan length | 468 mm |
| | 80 kV |
| CTDI _{vol} | 3.04 mGy |
| DLP | 149 mGy cm |

- Excellent image quality at low dose thanks to High Power 80
- 14 year old adolescent




3.0 mm MPRs

Pediatric Imaging



| | |
|---------------------|----------|
| Scan time | 0.9 s |
| Scan length | 162 mm |
| | 70 kV |
| CTDI _{vol} | 0.43 mGy |
| DLP | 9 mGy cm |

- CARE Child protocols with 70 kV for low dose in pediatrics
- 5 months old with collapsed lung



Tin Filter Topogram

Pediatric Imaging



Tin Filter Topogram

| | |
|---------------------|-----------|
| Scan time | 0.8 s |
| Scan length | 138 mm |
| | 70 kV |
| CTDI _{vol} | 0.52 mGy |
| DLP | 10 mGy cm |

- Fast acquisition for high quality pediatric scanning
- Follow-up scan after intervention

Acute Care and Trauma


Acute Care and Trauma

| | |
|---------------------|------------|
| Scan time | 5 s |
| Scan length | 335 mm |
| | 90 kV |
| CTDI _{vol} | 10.46 mGy |
| DLP | 301 mGy cm |

- Lung imaging with 90 kV in case of pulmonary embolism
- Great image quality with 1 mm MPRs
- Low contrast media application of 50 ml / 350 mg/ml



MIP



1 mm MPRs




1 mm MPRs




Acute Care and Trauma

| | |
|---------------------|------------|
| Scan time | 3 s |
| Scan length | 390 mm |
| | 120 kV |
| CTDI _{vol} | 8.38 mGy |
| DLP | 292 mGy cm |

- High-resolution lung imaging thanks to Stellar detector and high number of projections
- Fast scan speed to adapt to patient's breathhold capabilities



1 mm MPRs



MIP



syngo.CT DE Lung Analysis

Spectral Imaging with Dual Energy

| | |
|---------------------|------------|
| Scan time | 9.7 s |
| Scan length | 382 mm |
| | AuSn120 kV |
| CTDI _{vol} | 4.52 mGy |
| DLP | 153 mGy cm |

- TwinBeam Dual Energy
- Global decrease of lung perfusion after pulmonary embolism
- Zero-click spectral imaging reconstructions with Recon&GO

Acute Care and Trauma

| | |
|---------------------|-------------|
| Scan time | 11 s |
| Scan length | 1011 mm |
| | 110 kV |
| CTDI _{vol} | 10.2 mGy |
| DLP | 1072 mGy cm |

- Excellent imaging quality even in challenging positioning with both hands down



1 mm MPR



MIP



1 mm MPR




Cinematic VRT

Acute Care and Trauma

| | |
|---------------------|------------|
| Scan time | 7.3 s |
| Scan length | 167 mm |
| | 110 kV |
| CTDI _{vol} | 22.53 mGy |
| DLP | 424 mGy cm |

- Sharp visualization of fracture thanks to improved resolution with Stellar detector and UFC
- Orbita fracture evaluation of two years old child



0.6 mm MPR



2 mm MPR



Cinematic VRT


Acute Care and Trauma

| | |
|---------------------|------------|
| Scan time | 8 s |
| Scan length | 207 mm |
| | Sn110 kV |
| CTDI _{vol} | 10.3 mGy |
| DLP | 169 mGy cm |

- Tin Filter with 10 kV Steps for excellent image quality



1 mm MPRs



Cinematic VRT

Acute Care and Trauma

| | |
|-----------|------|
| Scan time | 15 s |
|-----------|------|

| | |
|-------------|--------|
| Scan length | 487 mm |
|-------------|--------|

| | |
|--|--------|
| | 120 kV |
|--|--------|

| | |
|---------------------|----------|
| CTDI _{vol} | 6.61 mGy |
|---------------------|----------|

| | |
|-----|------------|
| DLP | 346 mGy cm |
|-----|------------|


- Adding diagnostic information with iterative Metal Artifact Reduction (iMAR) also in acute care and trauma settings



With iMAR



Without iMAR




With iMAR
1 mm slice MPRs



Musculoskeletal Imaging

Musculoskeletal Imaging

| | |
|---------------------|------------|
| Scan time | 9.2 s |
| Scan length | 439 mm |
| | Sn130 kV |
| CTDI _{vol} | 8.33 mGy |
| DLP | 393 mGy cm |

- Combination of high kV and Tin Filter scanning for powerful metal artefact surpression





Spectral imaging with Dual Energy


| | |
|---------------------|----------------|
| Scan time | 8 s |
| Scan length | 294 mm |
| | 100 / Sn140 kV |
| CTDI _{vol} | 6.84 mGy |
| DLP | 250 mGy cm |

- TwinSpiral Dual Energy
- Excellent spectral separation thanks to Tin Filter enables imaging of the bone marrow
- Zero-click postprocessing with inline results in Recon&GO

Musculoskeletal Imaging

| | |
|---------------------------|------------|
| Scan time | 14 s |
| Scan length | 329 mm |
| | Sn140 kV |
| CTDI_{vol} | 15.75 mGy |
| DLP | 451 mGy cm |

- High kV Tin Filter scanning for efficient metal artefact suppression
- Obese patient



Cinematic VRT


Musculoskeletal Imaging

| | |
|---------------------|------------|
| Scan time | 7 s |
| Scan length | 316 mm |
| | Sn110 kV |
| CTDI _{vol} | 4.68 mGy |
| DLP | 128 mGy cm |

- High resolution for great visualization of fracture
- Tin Filter for low-dose trauma scan at great image quality



Cinematic VRT



1 mm MPRs


CT-guided Intervention Guide&GO

CT-guided Intervention Guide&GO

Native and contrast enhanced
CT for planning (each):

| | |
|---------------------|------------|
| Scan time | 6 s |
| Scan length | 284 mm |
| | 80 kV |
| CTDI _{vol} | 13 mGy |
| DLP | 367 mGy cm |

- CT-guided biopsy with Guide&GO
- Intuitive, tablet-based workflow from planning CT to control scan



CT-guided Intervention Guide&GO

i-Sequence Abdomen (2x)
for biopsy:

| | |
|---------------------|----------------|
| Scan time | 0.7 / 2 s |
| Scan length | 16 mm |
| | 120 kV each |
| CTDI _{vol} | 6.8 / 8.1 mGy |
| DLP | 10 / 12 mGy cm |

- CT-guided biopsy with Guide&GO
- Safeguard needle placement with FAST i-sequence for instant monitoring of the needle position



How to get there

How to get there



myExam Companion

SOMATOM go. platform starts the era of intelligent CT scanning with myExam Companion. With it, AI turns aggregated data into built-in expertise to automatically leverage the full potential of technologies – regardless of where, when, and by whom the results have been produced.

This allows for reliable and reproducible results from day one. myExam Companion means performing all advanced CT examinations as perfectly and quickly as if they were routine, allowing even unexperienced users to find the best combination of parameters for every individual patient and procedure. Personalized imaging for precise dose and contrast media optimization, reduced unwarranted variations, and always consistent results for enhanced diagnostic experience.




Inline/offline postprocessing

SOMATOM go.Top comes with two kinds of postprocessing tools: a zero-click “inline” reconstruction toolkit and another one for “offline” diagnosis.

The inline results of Recon&GO save time, reduce workflow steps, and deliver ready-to-read, standardized images. As a standard, Recon&GO includes anatomical ranges, table and bone removal, vascular CPR (Curved Planar Reconstruction), and multi recon (for automated multiple reconstructions in just one step). The High Performance Package adds spine ranges and rib ranges. Additionally it can be complemented with dedicated inline results for more advanced clinical tasks such as cardio, Dual Energy, neuro and pulmonary imaging.

For offline diagnosis, CT View&GO offers dedicated tools for smooth and efficient reading. Its standard version includes anatomical ranges, table and bone removal, vessel extension, spine ranges and endoscopic view among others. Furthermore, it can be extended with options in the field of cardio and Dual Energy.

Additionally, you can purchase *syngo*.CT Osteo, *syngo*.CT Neuro DSA and dedicated *syngo*.CT Dual Energy applications for different clinical questions.




Holistic Dual Energy solution

SOMATOM go.Top offers a holistic solution with two Dual Energy modes that makes DE routine ready for all patients and situations – completely neutral in dose and workflow.

TwinSpiral is supported by a new workflow concept of two scans integrated into one single acquisition. It offers the possibility to acquire two spiral data sets in sequence at different energies. Thanks to the spectral properties of the Tin Filter, TwinSpiral DE allows a better spectral separation for non-contrast examinations, whereas TwinBeam Dual Energy is especially useful for characterizing contrast media examinations since it acquires low and high-kV datasets in a single scan. By allowing you to characterize, highlight, and quantify different materials, Dual Energy gives you greater diagnostic confidence with virtually all patients.

And combined with dedicated Spectral Viewing packages it allows for comprehensive assessment. No matter if you would like to do your postprocessing directly at the AWP or rather have it sent automatically to PACS by Recon&GO – it offers a solution for all clinical workflows.



Athlon™ tube

High Power 70 allows you to scan at the highest tube current in its class: up to 825 mA at 70 kV. This achieves better iodine contrast for sharper images, even in small distal vessels. As a result, you can considerably reduce contrast media and thus scan more patients, deliver better patient care, and reduce examination costs.

Furthermore, the tube voltage is automatically tailored to each patient and clinical indication by CARE kV. Voltage levels can be adjusted at intervals of 10 kV for less dose and high contrast resolution and are aligned with respective tube currents. This keeps dose low, while image quality stays excellent.



Stellar detector

The Stellar detector reduces image noise in every scan, while the advanced iterative reconstruction SAFIRE² delivers excellent image quality at very low dose. This provides excellent and homogenous image quality, even in complex areas.

The Stellar detector's high-end technology includes fully integrated components and an advanced 3D antiscatter collimator. It keeps electronic noise low, increases dose efficiency, and improves spatial resolution. The smart configuration of the detector elements simplifies access, eases maintenance, and increases scanner uptime.



Tin Filter

Inherited from high-end dual source scanners, the Tin Filter (Sn) cuts out lower energies to reduce dose and optimizes contrast between soft tissue and air. This has direct benefits in lung and colon imaging, for example. Clinical experience also shows that Tin Filter technology reduces beam-hardening artifacts and improves image quality in bony structures, making it extremely useful in orthopedic examinations.



On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens Healthineers sales organization worldwide. Availability and packaging may vary by country and is subject to change without prior notice.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features which do not always have to be present in individual cases.

Siemens Healthineers reserves the right to modify the design, packaging, specifications, and options described herein without prior notice. Please contact your local Siemens Healthineers sales representative for the most current information.

Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

¹ Automated Landmarking and Parsing of Human Anatomy

² In clinical practice, the use of SAFIRE may reduce CT patient dose depending on the clinical task, patient size, anatomical location, and clinical practice. A consultation with a radiologist and a physicist should be made to determine the appropriate dose to obtain diagnostic image quality for the particular clinical task.

International version.

Not for distribution or use in the U.S.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone: +49 9131 84-0
siemens-healthineers.com