

Cutting-edge robotic imaging to drive minimally invasive procedures for multidisciplinary usage



| Speed | Precision | No repetition |
|---|---|--|
| Optimize diagnostic experience by saving time with software that supports you in planning your needle paths. <i>syngo</i> Needle Guidance enables navigation even without the help of an optical navigation system. ARTIS pheno also automatically moves into position for the planned screw path and guides you to the precise entry point with a laser cross. | Visualize all anatomical details and plan your procedure with a 4-second <i>syngo</i> DynaCT in the prone position. ARTIS pheno allows you to optimally treat nearly any patient, regardless of patient size. The power-assisted tabletop permits hassle-free repositioning with zero force, even in tilted orientations for patients up to 280 kg. | Improve clinical outcomes by offering immediate postoperative quality control with the 4-second <i>syngo</i> DynaCT and <i>syngo</i> Fusion Package. Maintain the highest infection control standards in the hybrid OR with the antimicrobial surfaces of ARTIS pheno and its floor-mounted design that ensures uninterrupted laminar airflow. |

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization world-wide.

Availability and packaging may vary by country and are subject to change without prior notice. Some/ALL of the features and products described herein may not be available in the United States.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Siemens reserves the right to modify the design, packaging, specifications, and options described herein without prior notice. Please contact your local Siemens sales representative for the most current information.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

The products/features and/or service offerings mentioned herein are not commercially available in all countries and/or for all modalities. If the services are not marketed in countries due to regulatory or other reasons, the service offering cannot be guaranteed. Please contact your local Siemens organization for further details.

The statements by Siemens' customers described herein are based on results that were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption) there can be no guarantee that other customers will achieve the same results.

¹ Courtesy: Jikei University School of Medicine, Tokyo, Japan

² Courtesy: Diakonie Klinikum Jung-Stilling, Siegen, Germany

³ Quote: Dr. Ohashi, Senior Spine Surgeon – Neurosurgery Dept., Jikei University School of Medicine, Tokyo, Japan

.....
Siemens Healthineers Headquarters
 Siemens Healthcare GmbH
 Henkestr. 127
 91052 Erlangen, Germany
 Phone: +49 9131 84-0
 siemens-healthineers.com

Published by Siemens Healthcare GmbH · Order No. A91AT-23205-60C1-7600 · Printed in Germany · 0000.x · ©Siemens Healthcare GmbH, 2018

siemens-healthineers.com



Deliver outcomes that matter to patients undergoing spinal fusion surgery, with ARTIS pheno

Plan and place pedicle screws precisely with *syngo* Needle Guidance



Speed. Precision. No Repetition.

Thanks to improving technologies, minimally invasive procedures are increasing in spine surgery, leading to reduced collateral tissue damage and scarring. This ultimately lowers healthcare costs and leads to better patient outcomes compared to open spine surgery.

To expand precision medicine, a reliable technology is crucial in spine surgery settings: image-guided or robotic-assisted systems are new partners that help make spine procedures **more safe, accurate, and efficient** by minimizing the risk of complications and increasing patient safety.

One step in this direction is the Siemens Healthineers-engineered ARTIS pheno: a cutting-edge robotic imaging system for individualized **preoperative planning, intraoperative guidance, and immediate postoperative quality control** – with *syngo* Needle Guidance you can even plan and place pedicle screws precisely without the help of a navigation system.



Better results in challenging spinal fusion cases

Case 1 Challenges:


- Precise screw placement
- Use of large instruments

Case 2 Challenges:


- Use of large instruments
- Obese patient (BMI 44)

ARTIS pheno offers:

- syngo DynaCT Body in 4 seconds
- Wide-space C-arm
- syngo Needle Guidance
- syngo DynaCT SMART for metal artefact reduction



Case 1¹




Case 2²




Preoperative Planning the screw path you need for surgery

Preoperative planning can be performed in the OR right before surgery with an intraoperative syngo DynaCT 3D scan in just 4 seconds. This saves time and the hassle of outside planning.

Plan the needle path using syngo Needle Guidance on a scan performed in the prone position (as opposed to a preoperative supine image), which produces an exact anatomical outline. An entry and a defined endpoint can be planned easily for multiple pedicle screw pathways. Screw lengths can be calculated and the OR team can plan the procedure accordingly.




syngo Needle Guidance: planning the procedure on an intraoperative syngo DynaCT scan¹



Intraoperative Deployment facilitated by syngo Needle Guidance

syngo Needle Guidance helps you define the pedicle screw paths. The laser cross then guides you to the desired entry point at the angle you envisioned. After defining the needle path using syngo Needle Guidance, the ARTIS pheno automatically moves in the planned direction, pinpointing the entrance point with a laser cross. The path is shown superimposed on the live fluoroscopy image regardless of the perspective you choose.

This can make a navigation system obsolete in 90% of cases³ – save for major deformities or complex procedures. In sum, the ARTIS pheno offers more space, increased speed and can help to reduce patient burden.



Inserting K-wires and screws under fluoroscopic guidance¹




Intraoperative Deployment using the wide-space C-arm

The wide-space C-arm with free space of 95.5 cm between the tube and detector allows working with large instruments during minimally invasive procedures. This makes patient positioning simpler.

Working with large instruments, especially for thoracic spine procedures, lateral positioning, or use of a microscope, can now be performed without pushing the C-arm out of the way or sacrificing workspace.




Operating with multiple products with the C-arm still in place: lateral a.p.¹



Intraoperative Deployment of implants using low dose fluoroscopy

It is most important to ensure immediate correct placement of the pedicle screws or other implants used. High-end fluoroscopy image quality enables precise positioning even at low dose.




2D fluoroscopy image¹



Postoperative Quality control with a 4-second syngo DynaCT

At the end of surgery, right before the patient is closed, a postoperative syngo DynaCT scan of 4 seconds can confirm precise positioning of the implants used. This makes revision surgeries to correct misplaced pedicle screws obsolete, lowers the burden on the patient, and reduces healthcare costs.



Postoperative syngo DynaCT¹

Postoperative syngo DynaCT²



“One huge benefit is the additional 13 cm of SID, which helps in performing a syngo DynaCT without colliding with the head clamp, or in performing thoracic spinal fusion, which requires large instruments.”

Prof. Yuichi Murayama, MD,
Jikei University School of Medicine, Tokyo, Japan



“Since we have ARTIS pheno in our OR, our colleagues really enjoy their work.”

Prof. Veit Braun, MD,
Diakonie Klinikum Jung-Stilling, Siegen, Germany