

## Study Protocol

# Peripheral angiography using CO<sub>2</sub>

Interventional Radiology

Artis zeego with variable isocenter is very helpful in CO<sub>2</sub> cases as it travels easily along the tilted table.

---

### Courtesy of

Ulf Teichgräber, MD,  
Renè Aschenbach, MD,  
Department of Diagnostic and  
Interventional Radiology,  
Jena University Hospital,  
Germany

### System & Software

Artis zeego with Q technology VD10  
syngo X Workplace VC10

---

## Case Description

### Patient history

74-year-old female patient. Occlusive peripheral artery disease. Patient did not want bypass surgery. Renal insufficiency Grade III (severely limiting the use of contrast media).

### Diagnosis

Fontaine IIb. Left superficial femoral artery is occluded in the proximal third and there is distal filling via collaterals.

### Treatment

Laser atherectomy followed by PTA with drug-eluting balloon. VIABAHN stent graft placed.

### General comments

CO<sub>2</sub> evenflow acquisition technique provides high image quality. High k-factor creates more homogeneous, bubble-free visualization of vessels.

## Peripheral angiography using CO<sub>2</sub>

Acquisition protocol	DSA CO <sub>2</sub> Evenflow
Frame rate	7.5 f/s (higher frame rate, low dose/f, high k-factor)
Length of sequence	20 s max
<b>Injection protocol</b>	
Catheter position	Distal abdominal artery
Contrast medium (CM)	CO <sub>2</sub> 1.3 bar
Dilution (CM/Saline):	No
Injection volume	80 mL for each angiography step/for intervention 20 mL CO <sub>2</sub>
Injection rate	Manual
Duration of injection	Manual
X-ray delay	Injection starts right after mask selection (~4th frame)
Power injector used	Optimed CO <sub>2</sub> system
<b>Image Postprocessing</b>	
Pixel Shift	Automatic
Min OPAC	Yes

## Clinical Images




Figure 1: Min OPAC – Iliac arteries




Figure 2: Min OPAC – Vessels below occlusion show collateral filling




Figure 3: Axial MPR – primary reconstruction




Figure 4: Axial MPR – secondary reconstruction

---

### Siemens Healthineers Headquarters

Siemens Healthcare GmbH

Henkestr. 127

91052 Erlangen, Germany

Phone: +49 9131 84-0

siemens-healthineers.com

*The statements by Siemens' customers presented here are based on results that were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption), there can be no guarantee that other customers will achieve the same results.*

*On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this case are available throughout the Siemens sales organization worldwide.*

*All rights reserved.*