

SOMATOM On.site

# Bringing critical care imaging to your patient

Focus: Radiology

[siemens-healthineers.com/somatom-onsite](https://siemens-healthineers.com/somatom-onsite)



reddot winner 2020  
best of the best

International version.  
Not for distribution or use in the U.S.

**SIEMENS**  
Healthineers



3.202-3.224  
3.225-3.24

Hill-Rom

# Providing CT scans of the head for ICU patients is crucial but often challenging

ICU (Intensive Care Unit) patients suffering from acute and critical head conditions face a high risk of sudden health status deterioration. When this happens, fast evaluation is essential for immediate treatment decision-making. CT scans are ideally suited for this task – but providing imaging for ICU patients is often cumbersome and not without risk.



## CT scans of the head are in high demand

Head scans are among the most common procedures for ICU patients and account for 30% of all CT scans.<sup>3</sup>



## ICU patient transport involves logistical and staffing challenges

Transporting patients from the ICU to the radiology department is a delicate and complex process that requires up to five staff members<sup>1</sup> to leave their other duties behind.




## ICU staff shortage limits patient care

Frequent bed closures affect 76% of German ICU floors, mostly due to staff shortage.<sup>4</sup>



## Transport may negatively affect the patient's health status

Providing CT imaging requires transporting ICU patients to the radiology department. However, for 1 in 5 patients, transport makes existing health conditions worse.<sup>2</sup>



## ICU patient scans complicate radiology department schedules

Patients with acute and critical conditions tie up CT scanners for at least twice as long as most other patients.<sup>5,6</sup> Transports are also frequently delayed,<sup>2</sup> making it necessary to reschedule other patients.





## Bringing critical care imaging to your patient

Our mobile head CT scanner SOMATOM On.site provides access to consistent and reliable SOMATOM image quality, right at your patient's bedside. Thanks to the latest advances in digital healthcare, the intelligent software myExam Companion together with our GO technologies guides users through the scanning process. With this workflow guidance the system is easy to use and it enables greater consistency in image acquisition across your team.

SOMATOM On.site will change the way you care for ICU patients suffering from acute and critical head conditions. Reducing transport to the radiology department may help prevent related complications and allows your staff to concentrate on what matters most: providing optimal care for every patient at the ward. With this easy access to CT imaging of the head directly at the point of care, you can get images to help diagnose and treat sudden changes in the patient's condition or arising complications faster. Additionally, offering CT imaging directly in the ICU allows the radiology department to maintain their scheduling.

**With SOMATOM On.site, you can transform care delivery.**



### **See for yourself**

Our workflow video demonstrates how point-of-care CT scanning works with SOMATOM On.site.

---

# Transform care delivery with SOMATOM On.site

## Reduce patient transports

SOMATOM On.site permits CT imaging of the head directly in the ICU ward, so patients need not encounter any transport-associated risks. Acquiring images directly at the patient's bedside may speed up diagnostics and following decisions and steps.

## Optimize staff deployment

SOMATOM On.site reduces the need for transporting patients to the radiology department for CT scans of the head. ICU staff and doctors can stay in the ICU. This means they can gain additional time for providing bedside patient care and focus more on their core tasks.


## Achieve reliable and consistent image quality at the point of care

SOMATOM On.site brings SOMATOM image quality to the patient's bedside. It is combining proven technologies from existing SOMATOM scanners with newly designed components.


Point-of-care imaging reduces the number of ICU patient scans that need to be performed in the radiology department, potentially resulting in more streamlined utilization of stationary systems and staff, along with more predictable radiology appointment schedules.




# Imaging at the patient's bedside



Non-contrast CT images of the head with ALPHA\* anatomy recognition technology




Delayed contrast media-enhanced examinations




All orientations and reformations PACS-ready

## Apply known tools and applications



Maximum intensity projection (MIP), e.g., for intervention planning



Cinematic Rendering to enhance education and communication

\* ALPHA = Anatomical Landmarking and Parsing of Human Anatomies

# Bring SOMATOM quality to the patient's bedside

## SOMATOM imaging performance

SOMATOM On.site delivers reliable and consistent image quality on a mobile CT scanner thanks to a combination of existing innovations from SOMATOM scanners (e.g., low-noise detectors) and newly designed components (e.g., the tube).



### Stellar detector

SOMATOM On.site is the first mobile CT system outfitted with our SOMATOM Stellar detector. Stellar detector's technology helps to keep electronic noise low, increases dose efficiency, and improves spatial resolution. It's, high channel density and integrated design of the detector results in the grey-white matter differentiation that is crucial for evaluating CT scans of the head.



### Dedicated X-ray tube

Our X-ray tube was designed specifically for neuro CT imaging. It features the resistance to vibrations and shock that is needed in mobile CT scanners. It is built to reduce maintenance downtime and improve performance in demanding environments. The cooling system is a special design combining water and air.


## State-of-the-art image reconstruction algorithms

SOMATOM On.site incorporates iterative reconstruction algorithms to help you get the most out of every scan.



### Metal artifact reduction with iMAR

iMAR (iterative metal artifact reduction) allows for reconstructed artifact reduction without losing valuable information. iMAR is designed to yield images with a reduced level of metal artifacts compared to conventional reconstruction if the underlying CT data is distorted by the presence of metal in the scanned object.



### Noise removal with SAFIRE

This iterative reconstruction algorithm delivers natural looking images at the lowest possible dose. Multiple iterations are performed when processing the raw data. During these iterations the noise is reduced without compromising image quality.

## Telescopic gantry for hassle-free acquisitions

SOMATOM On.site features a specially designed telescopic gantry with two independent gantry compartments. The rear compartment houses the tube and the detector and moves backwards during scan acquisition. The front part stays in place during scan acquisition, so the many lines and drains typically present in ICU patients, are not affected during scanning with SOMATOM On.site.

Thanks to this technology, the systems' trolley position remains fixed throughout the scan, eliminating common problems associated with scanners that may move during image acquisition. Taking scanner, table, and patient movement out of the equation means that you do not have to make any compromises when it comes to image quality.



Scan start position



Scan end position

## Technical specifications

### Key data

Slices	32
Z-coverage	24 mm (32 x 0.75 mm)
Detector material	Stellar detector UFC (Ultra Fast Ceramics)
Rotation time	1.0 s
Max. scan speed	36 mm/s

kV settings	80kV; 120 kV
Current	up to 25 mA
Scanner dimensions (height / width / depth)	1.55 / 1.60 / 0.74 m
Iterative reconstruction	SAFIRE, iMAR
Scanner weight	~ 998 kg



## Patient positioning concept

SOMATOM On.site features a telescopic gantry and stationary system trolley that remains retracted until you are ready to start scanning. This gives you ample space for patient positioning with tubes and lines in place. The scanner comes with convenient integrated patient support accessories. The integrated body support (shoulder board for neuro imaging) helps bridging the gap between

the bed and scanner. Sliding the patient's upper body onto the shoulder board permits convenient patient positioning in preparation for scans. Furthermore, the integrated patient scan support (head holder for neuro imaging) of SOMATOM On.site helps for correct isocentric positioning of the patient's head inside the gantry bore.



## Radiation safety concept

Traditional stationary CT scanners require users to return to the control room during scans, because a safety distance is necessary for radiation protection. SOMATOM On.site allows you to stay close to patients and remain at their side. The self-shielded scanner design, along with the detachable radiation shields covering the front and back of the gantry bore, reduce scatter radiation and provide protection for staff and neighboring patients.



## Fast and intuitive scan preparation

SOMATOM On.site leverages key digital healthcare innovations to unleash the full potential of modern imaging technology – independent of the operator’s experience. The scanner offers an intuitive user interface on the Touch UI and, thanks to myExam Companion, guides you through the complete imaging examination.

In addition SOMATOM On.site features workflow guidance provided by our proven GO technologies. Scan&GO enables techs to control the entire scan process via the integrated Touch UI and stay right next to the patient throughout the image acquisition process. After the scan the users have the possibility to check the images (Check&GO) before sending them to PACS. A dedicated Recon&GO functionality performs automated postprocessing and uploads the images to PACS without any user interaction. As a result, radiologists receive comprehensive results and consistent image quality. Together, these features make scanning intuitive and easy.

## Simple steps suitable for cross-trained staff



Select your patient via the intuitive user interface of SOMATOM On.site.



Choose your preferred scan protocol for image acquisition.



Scan&GO allows you to supervise the scan right at your patient’s bedside.



Check&GO gives you the chance to check images before sending them to PACS.



Recon&GO handles all postprocessing for you and automatically uploads the images to the PACS.

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens Healthineers sales organization worldwide.

Availability and packaging may vary by country and is subject to change without prior notice. Some or all of the features and products described herein may not be available in the United States.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Siemens Healthineers reserves the right to modify the design, packaging, specifications, and options described herein without prior notice. Please contact your local Siemens Healthineers sales representative for the most current information.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

1 Barnes E. *Portable CT brings lower costs, reduce risk to ICU.* AuntMinnie.com; 2008.

2 *Anticipating unique disruptions. PwC's 21st CEO Survey: key findings from the healthcare industry [homepage on the internet].* PwC Turkey; 2018 [cited 2019 Aug 6]. Available from: <https://www.pwc.com.tr/en/yaayinlar/arastirmalar/ceo-survey/21-ceo-arastirmasi/saglik-sektoru-sonuclari.html>.

3 Data on file.

4 Karagiannidis C, Kluge S, Riessen R, Krakau M, Bein T, Janssens U. *Auswirkungen des Pflegepersonalmangels auf die intensivmedizinische Versorgungskapazität in Deutschland.* *Med Klin Intensivmed Notfmed.* 2018 Jul 9;114(4):327–33. Available from: <https://www.springermedizin.de/auswirkungen-des-pflegepersonalmangels-auf-die-intensivmedizinis/15927204>.

5 Agrawal, D et al. (2016). *Bedside computed tomography in traumatic brain injury: Experience of 10,000 consecutive cases in neurosurgery at a level 1 trauma center in India.* *Neurol India.* 2016 Jan–Feb; 64(1):62–365. Available from: <https://www.ncbi.nlm.nih.gov/pubmed/26754994>.

6 Masaryk et al. (2008). *The economic and clinical benefits of portable head/neck CT imaging in the intensive care unit.* *Radiol Manage.* 2008 March–April; 30(2):50–54. Available from: <https://www.ncbi.nlm.nih.gov/pubmed/18431942>.

---

### **Siemens Healthineers Headquarters**

Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen, Germany  
Phone: +49 9131 84-0  
[siemens-healthineers.com](http://siemens-healthineers.com)