



Analytical and Clinical Performance of the ADVIA Centaur Anti-thyroglobulin (aTgII) Assay*

Arcuri M, Low O, Krumm C, Ray S, Dauscher C.
Siemens Healthcare Diagnostics Inc., Tarrytown, NY, U.S.

[siemens-healthineers.com](https://www.siemens-healthineers.com)

Clinical
Brief

* Not available for sale in the U.S. Product availability varies by country and is subject to regulatory requirements.

Introduction

Thyroglobulin (Tg) is a large, heterogeneous glycoprotein found in the follicular cells of the thyroid. Thyroglobulin plays an important role in the biosynthesis of thyroid hormones triiodothyronine (T3) and thyroxine (T4). Within the thyroid gland, the enzyme thyroid peroxidase catalyzes the iodination of thyroglobulin's tyrosyl groups. The iodinated thyroglobulin is stored in the colloid of the follicle and serves as a storage reservoir for T3 and T4. When the thyroid gland is stimulated, thyroglobulin is degraded, and energy-regulating hormones T3 and T4 are released into the bloodstream.¹⁻³

Production of autoantibodies against thyroglobulin (anti-Tg) is a result of the immune system targeting components of the thyroid gland. The measurement of autoantibodies against thyroglobulin (anti-Tg) is useful in differentiating patients with autoimmune thyroid disease.⁴⁻¹³ In patients with thyroid carcinoma, measurement of Tg antigen must consider the likelihood of the presence of significant levels of anti-Tg since measurement and detection of Tg antigen may be influenced by the presence of anti-Tg antibodies.¹⁴⁻¹⁵

Assay Principle

The ADVIA Centaur® Anti-thyroglobulin II (aTgII) Assay* is a fully automated, one-step, analyte-bridging immunoassay using acridinium ester chemiluminescent technology for the detection of thyroglobulin autoantibodies. This assay uses human thyroglobulin in both the Lite Reagent and the Solid Phase. A direct relationship exists between the amount of anti-Tg present in the patient sample and the number of light units (RLUs) detected by the system.




Figure 1. ADVIA Centaur aTgII Assay format.

Table 1. Assay attribute summary.

Sample Volume	25 µL
Time to Result	18 minutes
Specimen Types	Serum, plasma (lithium heparin and EDTA)
Measuring Interval	1.8–1000 IU/mL
Limit of Detection	1.3 IU/mL
Limit of Quantitation	1.8 IU/mL
Hook Effect	No hook effect up to 50,000 IU/mL

*Not available for sale in the U.S. Product availability varies by country and is subject to regulatory requirements.

Experimental Design

Analytical performance was studied in terms of precision, detection capability (LoB, LoD, and LoQ), interference, linearity, and hook effect. Clinical performance was established using a reference interval study and clinical concordance against a commercially available assay.

Limit of quantitation (LoQ)

- CLSI protocol EP17-A2¹⁶
- Nine serum pools, two runs/day for 5 days
- One reagent lot
- Within-laboratory precision calculated by analysis of variance (ANOVA)
- LoQ determined as the concentration at which the within-laboratory CV is $\leq 20\%$ CV

Linearity

- CLSI protocol EP06-A¹⁷
- 11 samples, tested $n = 3$
- One reagent lot
- Analysis by percent bias from weighted linear regression

Precision

- CLSI protocol EP05-A3¹⁸
- Six specimens, two runs/day for 20 days
 - Six patient serum pools
- One reagent lot
- Data calculated by analysis of variance (ANOVA)
- Data reported as the upper 95% confidence limit

Interference

- CLSI protocol EP07-ed3¹⁹
- Three analyte test levels tested: 5, 50, and 500 IU/mL, tested $n = 5$
- Analysis by percent bias

Reference interval

- 198 serum specimens from apparently healthy males <30 years
- Tested in singlicate on one ADVIA Centaur XP system and one Beckman ACCESS 2 system
- One reagent lot on each system
- Analysis by receiver operating characteristic (ROC) curve

Assay clinical cutoff vs. predicate

- CLSI protocol EP09-ed3²⁰
- 293 patients (Hashimoto's disease, Graves' disease, anemia, thyroid disease, non-thyroid autoimmune disease, and healthy patients)
- Tested in singlicate on one ADVIA Centaur XP system and one Beckman ACCESS 2 system (TgAb assay)
- One reagent lot on each system
- Analysis by receiver operating characteristic (ROC) curve

Results

Limit of quantitation

A limit of quantitation of 1.8 IU/mL was achieved for the assay.




Figure 2. LoQ determination.

Assay linearity

Linearity of the assay was achieved within the measuring interval of 1.8–1000 IU/mL. The percent bias from the weighted linear fit was <10% across the assay range.




Figure 3. Assay linearity.

Precision

The precision of the assay met our design requirements. The 95% upper confidence limit for repeatability ranged from 1.8 to 6.0% CV. The 95% upper confidence limit for within-lab precision ranged from 4.2 to 7.7% CV.

Table 3. Precision.

Specimen	Dose (IU/mL)	Repeatability		Within Lab	
		SD	%CV	SD	%CV
1	6.73	0.4	5.9	0.5	7.3
2	13.0	0.8	6.0	1.0	7.7
3	18.3	0.7	4.0	1.0	5.6
4	50.4	1.6	3.1	2.4	4.8
5	496	8.9	1.8	20.7	4.2
6	872	27.9	3.2	49.3	5.7

Interferents

The ADVIA Centaur aTgII Assay was tested for interference from the substances listed in Table 4. No significant interference (<6%) was observed for these substances.

Table 4. Potential interferents.

Acetaminophen	Insulin
Aspirin	Iodine
Bilirubin (conjugated)	Lipemia using INTRALIPID
Bilirubin (unconjugated)	Low protein
Biotin at 3500 ng/mL	L-thyroxine (T4)
Hemoglobin	Methimazole
High protein	RF
Ibuprofen	T3 antibodies

High-dose hook effect

High anti-Tg concentrations can cause a paradoxical decrease in the RLU (high-dose hook effect). In this assay, patient samples with anti-Tg concentrations above the measuring interval and as high as 50,000 IU/mL will report >1000 IU/mL.

Assay reference interval




Figure 4. High-dose hook effect.

A reference interval for the presence of antithyroglobulin antibodies was established using 198 apparently healthy males <30 years of age who had no personal or family history of thyroid disease and no history of autoimmune disease. The 95% nonparametric upper reference limit was calculated to be <1.3 IU/mL. A cutoff value for autoimmune thyroid disease was determined using a receiver operating characteristic (ROC) curve to optimize clinical agreement with the Beckman anti-Tg II assay.

Clinical cutoff vs. predicate




Figure 5. Nonparametric results (reference level 95%) for male patients <30 years old.

A total of 293 individuals were included in the study. With a cutoff value of 4.5 IU/mL, the positive percent agreement was 98.6% (71/72), with a 95% confidence interval (CI) of 92.5–100.0%. The negative percent agreement was 94.6% (209/221), with a 95% CI of 90.7–97.2%. The overall percent agreement was 95.6% (280/293), with a 95% CI of 92.5–97.6%.

Table 5. Method comparison.

ADVIA Centaur aTgII Assay	Beckman ACCESS 2 Tg-ab assay		
	Positive (≥4.0 IU/mL)	Negative (<4.0 IU/mL)	Total
Positive (≥4.5 IU/mL)	71	12	83
Negative (<4.5 IU/mL)	1	209	210
Total	72	221	293

Conclusion

The ADVIA Centaur aTgII Assay combines strong analytical performance with the workflow advantages of the automated ADVIA Centaur Immunoassay Systems.

- Fully automated assay on the ADVIA Centaur system with results in 18 minutes for serum and plasma samples
- Cutoff for autoimmune disease of 4.5 IU/mL
- LoQ of 1.8 IU/mL
- Repeatability $\leq 6.0\%$; within-lab precision $\leq 7.7\%$
- No interference to 3500 ng/mL biotin
- High positive (98.6%), negative (94.6%), and overall 95.6% agreement with the Beckman ACCESS Thyroglobulin Antibody II assay

At Siemens Healthineers, our purpose is to enable healthcare providers to increase value by empowering them on their journey toward expanding precision medicine, transforming care delivery, and improving patient experience, all made possible by digitalizing healthcare.

An estimated 5 million patients globally benefit every day from our innovative technologies and services in the areas of diagnostic and therapeutic imaging, laboratory diagnostics, and molecular medicine, as well as digital health and enterprise services.

We are a leading medical technology company with over 120 years of experience and 18,000 patents globally. Through the dedication of more than 50,000 colleagues in 75 countries, we will continue to innovate and shape the future of healthcare.

ADVIA Centaur and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc., or its affiliates. All other trademarks and brands are the property of their respective owners.

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

References:

1. Feldt-Rasmussen U. Analytical and clinical performance goals for testing autoantibodies to thyroperoxidase, thyroglobulin, and thyrotropin receptor. *Clin Chem.* 1996;42(1):160–3.
2. Torr ns JI, Burch HB. Serum thyroglobulin measurement: utility in clinical practice. *The Endocrinologist.* 1996;6(2):125-44.
3. Lucas M, Fern ndez-Ulloa M. Thyroid. In: Kaplan LA, Pesce AJ, eds. *Clinical chemistry: theory, analysis, correlation.* 3rd ed. St. Louis: CV Mosby; 1996:868-72.
4. Rose NR, Burek CL. Autoantibodies to thyroglobulin in health and disease. *Appl Biochem Biotechnol.* 2000; 83(1–3):245-51.
5. Kasagi K, Kousaka T, Higuchi K, et al. Clinical significance of measurements of antithyroid antibodies in the diagnosis of Hashimoto's thyroiditis: comparison with histological findings. *Thyroid.* 1996;6(5):445-50.
6. Rosenbaum D, Davies TF. The clinical use of thyroid autoantibodies. *The Endocrinologist.* 1992 Jan;2(1):55-62.
7. Burek CL, Rose NR. Thyroglobulin autoantibodies. In: Peter JB, Shoenfeld Y, ed. *Autoantibodies.* Amsterdam: Elsevier Science B.V.; 1996:810-15.
8. RA, Gilbert FI Jr, Miyamoto LA, et al. The superiority of antimicrosomal over antithyroglobulin antibodies for detecting Hashimoto's thyroiditis. *Arch Intern Med.* 1993 Apr 12;153(7):862-5.
9. Ruf J, Feldt-Rasmussen U, Hegeds L, et al. Bispecific thyroglobulin and thyroperoxidase autoantibodies patients with various thyroid and autoimmune diseases. *J Clin Endocrinol Metab.* 1994 Nov;79(5):1404-9.
10. Scherbaum WA. On the clinical importance of thyroid microsomal and thyroglobulin antibody determination. *Acta Endocrinol Suppl (Copenh).* 1987;281:325-9.
11. Walker DJ, Griffiths M, Griffiths ID. Occurrence of autoimmune diseases and autoantibodies in multi-case rheumatoid arthritis families. *Ann Rheum Dis.* 1986 Apr;45(4):323-6.
12. Ericsson UB, Christensen SB, Thorell JI. A high prevalence of thyroglobulin autoantibodies in adults with and without thyroid disease as measured with a sensitive solid-phase immunosorbent radioassay. *Clin Immunol Immunopathol.* 1985 Nov;37(2):154-62.
13. Weetman AP, McGregor AM. Autoimmune thyroid disease: further developments in our understanding. *Endocr Rev.* 1994 Dec;15(6):788-830.
14. Feldt-Rasmussen U, Rasmusen K. Serum thyroglobulin (Tg) in presence of thyroglobulin autoantibodies (TgAb). Clinical and methodological relevance of the interaction between Tg and TgAb in vitro and in vivo. *J Endocrinol Invest.* 1985;8:571-6.
15. Schaadt B, Feldt-Rasmussen U, Rasmusson B, et al. Assessment of the influence of thyroglobulin (Tg) autoantibodies and other interfering factors on the use of serum Tg as tumor marker in differentiated thyroid carcinoma. *Thyroid.* 1995;5(3):165-70.
16. Clinical and Laboratory Standards Institute. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures, approved guideline. EP17-A2. Wayne PA: CLSI; 2012.
17. Clinical and Laboratory Standards Institute. Evaluation of the linearity of quantitative measurement procedures: a statistical approach; approved guideline. EP06-A. Wayne, PA: CLSI; 2003
18. Clinical and Laboratory Standards Institute. Evaluation of precision of quantitative measurement procedures; approved guideline—third edition. EP05-A3. Wayne, PA: CLSI; 2014.
19. Clinical and Laboratory Standards Institute. Interference testing in clinical chemistry; approved guideline—third edition. EP07-ed3. Wayne, PA: CLSI; 2018.
20. Clinical and Laboratory Standards Institute. Measurement procedure comparison and bias estimation using patient samples; approved guideline—third edition. EP09c-ed3. Wayne, PA: CLSI; 2018.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone: +49 9131 84-0
siemens-healthineers.com

Published by

Siemens Healthcare Diagnostics Inc.
Laboratory Diagnostics
511 Benedict Avenue
Tarrytown, NY 10591-5005
USA
Phone: +1 914-631-8000