

White Paper

Guideline for integrated platelet aggregometry

For hemostasis systems from Siemens Healthineers

siemens-healthineers.com/hemostasis



Platelets: Crucial for primary hemostasis

In physiological hemostasis, the accumulation of platelets on the damaged vessel helps to stop bleeding by forming a platelet plug. There exist two types of platelet abnormalities that promote bleeding: qualitative and quantitative platelet disorders.

Platelets are formed by megakaryocytes in bone marrow and are essential for primary hemostasis. In the secondary phase of hemostasis, platelets act as a reactive surface for clotting factors and increase their activity, by which a stable fibrin clot is formed.

Platelets are small discoid cells without a nucleus. Although they do not contain DNA, platelets can synthesize various proteins using their RNA. In addition, they contain granules loaded with various metabolites, signaling molecules, and proteins, which are released when the platelets are activated.

Primary hemostasis is fundamentally important in the hospital. About 7% of acutely hospitalized patients show defects in primary hemostasis, which significantly increases the risk of bleeding.

Platelets: Activation, adhesion, and aggregation

When a vascular injury occurs, platelets bind via von Willebrand factor (VWF) to the collagen fibers of the now-accessible subendothelial matrix.

Platelets' main receptor for VWF is the glycoprotein (GP) Ib-IX-V complex. Under high shear stress, VWF first binds to collagen via the A3 domain and, in the next step, via the accessible A1 domain to the thrombocytic GPIb-IX-V-receptor (platelet adhesion).


After platelets have adhered to the injured tissue, they are activated by binding collagen to the platelet GPVI receptor.

This leads to a shape change and increased activation of the adherent platelets, with subsequent secretion of granular substances and consecutive aggregate formation.

The released granule ingredients (e.g., ADP, serotonin) increase platelet activation. This leads to recruitment of platelets from circulation, which in turn causes a further buildup of the platelet aggregate.

This is followed by the release of arachidonic acid and the formation of thromboxane A₂ (TxA₂), which enhances activation by binding to a specific thromboxane receptor. TxA₂ also has a vasoconstrictor effect. Within the activation the GPIIb-IIIa receptor changes its conformation to an activated form, so that it is able to bind fibrinogen.

Activated platelets are primarily linked by fibrinogen bridges through the GPIIb-IIIa receptor, VWF, and the GPIb-IX-V receptors. At the end of secondary hemostasis, the platelet aggregate is stabilized with cross-linked fibrin by factor XIII.⁷



Thrombocytopathy: Dysfunction of the platelets

Thrombocytopathy is a platelet dysfunction that can occur on several levels. These disorders are numerous and range from reduced secretion/synthesis of granules to defects in surface receptors. Since many signaling pathways that follow activation of the numerous

surface receptors are not yet or only incompletely known, clear classification of many thrombocytopathies is often difficult.

The most frequent thrombocytopathy is drug-induced platelet dysfunction.^{8,9}



Inherited thrombocytopathy

- Much less frequent
- Can lead to severe bleeding
- Often not diagnosed
- Examples: Glanzmann thrombasthenia, Bernard-Soulier syndrome, delta (δ)-storage pool disease (Hermansky-Pudlak syndrome), α -storage pool disease (gray platelet syndrome)



Acquired thrombocytopathy

- Frequent
- Drug-induced, e.g., by nonsteroidal anti-inflammatory drugs (NSAIDs), acetylsalicylic acid, and antiplatelet drugs
- Diseases, e.g., liver cirrhosis, kidney disease, myeloproliferative syndrome (MPS), multiple myeloma

Classification of platelet dysfunction⁷

Disease	Molecular Pathogens	Aggregometry
Receptor deficiencies		
Glanzmann thrombasthenia	Lack of expression or dysfunction of the GPIIb-IIIa receptor	Absent agonist-induced aggregation except ristocetin
Bernard-Soulier syndrome	Decreased expression or dysfunction of the GPIb-IX-V receptor	Absent or markedly reduced ristocetin-induced aggregation
von Willebrand syndrome Type 2B	Increased affinity of the GPIb-IX-V receptor to von Willebrand factor	Increased agglutination with low concentration of ristocetin
Defect of collagen receptor	Lack of expression or function of the GPIa-IIa or GPVI receptor	Absent or limited collagen-induced aggregation
Granular defects		
α-granular defects (e.g., gray platelet syndrome)	Absence of α-granules, disrupted release of adhesion proteins (e.g., fibrinogen, VWF) and membrane proteins (e.g., P-selectin, GPVI)	Reduced ristocetin- and collagen-induced aggregation
δ-granular defects (e.g., Hermansky-Pudlak syndrome)	Absence of δ-granules	Decreased ADP- and epinephrine-induced aggregation, absent or delayed second phase, disaggregation; reduced arachidonic acid- and collagen-induced aggregation
Signal transduction defects		
Aspirin-like defect	Disruption of the cyclooxygenase pathway (disorders of COX-1 enzyme or thromboxane synthetase)	Absent arachidonic acid-induced aggregation Failure of second phase of epinephrine-induced aggregation

Influences of drugs on platelet function

Most common are platelet dysfunctions as a result of taking medication, e.g., platelet aggregation inhibitors such as acetylsalicylic acid, NSAIDs, P2Y₁₂ receptor antagonists (clopidogrel, prasugrel, ticagrelor, cangrelor), or GPIIb-IIIa antagonists (abciximab, eptifibatide, tirofiban).⁶

Acetylsalicylic acid-induced platelet inhibition

Acetylsalicylic acid is the most often-used substance to inhibit platelet aggregation. It is established for the long-term prophylaxis of cardiovascular events.

It selectively and irreversibly inhibits cyclooxygenase 1 (COX-1) and thus the formation of thromboxane from arachidonic acid. Activation of thrombocytes via this pathway is thereby prevented.


P2Y₁₂ receptor antagonist-induced platelet inhibition

ADP receptor antagonists inhibit the P2Y₁₂ receptor (reversibly or irreversibly, depending on the mechanism of action).

Because of this inhibition, the intracellular cAMP cycle can no longer be influenced via the ADP-dependent activation of the P2Y₁₂ receptor.

GPIIb-IIIa antagonist-induced platelet inhibition


GPIIb-IIIa receptor antagonists are a category of drugs that bind to the activated GPIIb-IIIa receptor and, as a result, inhibit platelet aggregation. They are frequently used in percutaneous coronary interventions (intravenous only).



Diagnostic strategy for coagulation disorders

The recommended strategy is to support the diagnosis of inherited platelet disorders. The diagnostic steps depend on the patient's own and family history.⁴

AWMF-recommended diagnostic strategy



AWMF: Arbeitsgemeinschaft der Wissenschaftlichen Medizinischen Fachgesellschaften

ISTH: International Society on Thrombosis and Haemostasis

BAT: Bleeding Assessment Tool

MPV: Mean platelet volume

IPF: Immature platelet fraction (immature platelets)

ITP: Immune thrombocytopenia

TTP: Thrombotic thrombocytopenic purpura

HUS: Hemolytic-uremic syndrome

NGS: Next-generation sequencing

Pre-analytics: Sample collection and preparation

Pre-analytics plays a particularly important role in platelet aggregation testing.


Venous blood draw	Sample handling issues	Interference 1	Interference 2
• Draw blood carefully and gently to avoid contamination with tissue factor and hemolysis. Avoid activation of the platelets by shear stress. • Minimize use of tourniquets. • Preferred needle size is 21 gauge or larger. • Use only plastic or siliconized sample tubes with a citrate anticoagulant (3.2% or 105–109 mMol). • Do not use vacuum tubes if possible. • Fill the sample tube completely. • After collection, mix sample tube gently.	• Poorly mixed samples: Can lead to activation of coagulation factors and formation of small amounts of thrombin that can preactivate platelets. • Too-strongly shaken samples: Due to increased mechanical stress, preactivation of the platelets may occur. • Too-low transport temperature: Can lead to cryoactivation of the platelets. • Too-high transport temperature: Can lead to reduced platelet and sample stability.	• Hemolytic samples: Release of ADP, calcium, and hemoglobin leads to preactivation of platelets and influences optical measurement. • Lipemic samples: Increased optical density causes interference with optical measurement of platelet aggregation. • Icteric samples: Do not interfere with optical measurement of platelet aggregation. • Interferences can vary from system to system.	• Increased consumption of coffee and black tea as well as smoking influence platelet aggregation. • Stress leads to the release of epinephrine, which can cause preactivation of platelets. • Certain foods and food supplements can have an inhibitory or activating effect on platelet aggregation.⁴ • Antibiotic, psychotropic, and chemotherapeutic drugs influence platelet aggregation.^{8,9}

Platelet aggregometry: Sample preparation


For proper performance of platelet aggregation testing according to Born, follow these instructions for sample handling:

- To exclude thrombocytopenia, platelets should be determined in citrated whole blood, not EDTA whole blood. Use the THROMBOEXACT tube to most accurately determine the platelet count.
- Platelets are enriched during preparation of platelet-rich plasma (PRP) from citrated whole blood.
- An increased platelet count in the PRP indicates that the preparation was successful. A decreased platelet count in the PRP may indicate a preactivated sample.
- Allow the PRP to rest for 15–45 minutes at room temperature before starting the analysis. Gently swirl the PRP and place it on the system. Do not vortex the sample.
- The platelet count should be included in the diagnosis. If the platelet count decreases <100 Gpt/L, there is an increased influence on aggregation.
- Measurement of platelet aggregation should be completed within 2–4 hours after collection.

Platelet-rich plasma collection (PRP)



Platelet-poor plasma collection (PPP)



Platelet aggregometry: Light transmission aggregometry (LTA) according to Born


Light transmission aggregometry (LTA) according to Born⁷ is considered the gold standard of platelet function diagnostics. It can be used to determine and characterize many functional properties of platelets. Determination, interpretation of the findings, and assessment to define a clinical picture are very challenging.

In LTA, the increase in light transmission is determined over time after the addition of a specific agonist to the initially turbid PRP.

With the aid of the platelet-poor plasma (PPP), light transmission can be converted to % of aggregation. To obtain a standardized aggregation result, it is important to constantly stir the PRP to simulate shear stress. The temperature during aggregation measurement must be 37°C.

Integrated platelet aggregation on a routine hemostasis analyzer helps make platelet aggregation more manageable, reliable, and accessible, with easy-to-use, standardized testing.

Available parameters to characterize the aggregation curve



Up to nine different evaluation parameters are available on CN-6000 and CN-3000 Systems* as well as CS-2500 and CS-5100 Systems.

*The products/features/applications mentioned here are not commercially available in all countries and are subject to local regulations. Their future availability cannot be guaranteed. Not available for sale in the U.S.

Platelet aggregation agonists⁷

ADP	ADP activates the P2Y₁₂ receptor on the platelet surface and leads to inhibition of intracellular cAMP production, which results in a release of calcium. Furthermore, ADP activates the P2Y₁ receptor, initiating a platelet shape-change and a transient, rapidly reversible aggregation of the platelets. Activation of both receptors leads to a complete aggregation and release of contents from the platelet granules. Using a low ADP concentration, the aggregation may be biphasic or partially reversible.
Epinephrine	Epinephrine activates the platelets via the α_2-adrenoreceptor, and a biphasic aggregation curve occurs. Activation of the platelets by epinephrine leads to a platelet shape-change and a first wave of aggregation. As a result, fibrinogen receptors are exposed, and the activity of adenylatecyclases increases. A second wave of aggregation, quite different from the first, is caused by the release of ADP from the α-granules and the synthesis of thromboxane A₂.
Arachidonic acid	In the platelet, cyclooxygenase (COX-1) forms thromboxane A₂ from arachidonic acid. The TxA₂ leads to platelet activation over the thromboxane receptor. Platelet aggregation results in a simple aggregation wave.
Collagen	Collagen binds to the GPIIb-IIIa and GPVI receptors on the platelet surface. The GPVI receptor is partly responsible for platelet activation, which leads to the release of ADP, thromboxane A₂, and intracellular calcium. The aggregation curve typically shows a lag phase as a response to collagen, which is influenced by the collagen concentration. Platelet aggregation results in a simple aggregation wave.
Ristocetin	Ristocetin enables the interaction between von Willebrand factor and the glycoprotein GPIIb-V-IX receptor. By using low-dose ristocetin concentrations during aggregation, VWF type 2B can be detected.
TRAP	TRAP activates the thrombin receptor PAR-1 independently of receptor cleavage and mimics the effect of thrombin. TRAP-induced aggregation is not or only slightly inhibited by acetylsalicylic acid or P2Y₁₂ antagonists. TRAP can therefore be used to check whether platelets respond at all to a strong stimulus. PAR-1 antagonists show a concentration-dependent inhibition of aggregation with TRAP.

Aggregation: Integrated, automated, and standardized

On Siemens Healthineers mid- and high-volume hemostasis systems, many manual steps can be automated with the help of integrated platelet aggregation.* A high level of automation ensures reliable, reproducible, and safe processing of samples.

The following systems offer this integrated method:*

CN-6000 and CN-3000 Systems*

CS-5100 and CS-2500 Systems

To set up the method on your system, please contact your local Siemens Healthineers representative.

The example curves and explanations on the following pages assist in interpretation of test results.



CN-6000 and CN-3000 Systems*



CS-5100 System



CS-2500 System

**The products/features/applications mentioned here are not commercially available in all countries and are subject to local regulations. Their future availability cannot be guaranteed. Not available for sale in the U.S.*


Platelet aggregation integrated on the coagulation analyzer (1/2)*

Expected aggregation curves in normal patients

CN-6000 and CN-3000 Systems*

CS-2500 and CS-5100 Systems

Aggregation (%)



**The products/features/applications mentioned here are not commercially available in all countries and are subject to local regulations. Their future availability cannot be guaranteed. Not available for sale in the U.S.*


Platelet aggregation integrated on the coagulation analyzer (2/2)*

Expected aggregation curves in patients taking acetylsalicylic acid

CN-6000 and CN-3000 Systems*

CS-2500 and CS-5100 Systems

Aggregation (%)



*The products/features/applications mentioned here are not commercially available in all countries and are subject to local regulations. Their future availability cannot be guaranteed. Not available for sale in the U.S.

Notice

Results of light transmission aggregation (LTA), in particular if platelet counts of the PRP are lower or higher than recommended, should always be evaluated in conjunction with clinical history, clinical presentation, and other laboratory findings (such as bleeding time, complete blood count [CBC], and coagulation parameters).

In cases where aggregation results from tests with standard agonist concentrations do not agree with the clinical assessment, additional aggregation tests with lower agonist concentration should be performed. When using lower agonist concentrations, it is recommended to measure a healthy patient in parallel to evaluate the patient result. Moreover, further platelet function tests, e.g., ATP release, serotonin release, or flow cytometric assays, might be necessary.

Acknowledgements

We would like to thank all the experts who contributed to this guide.

Special thanks to:

Dr. med. Frauke Bergmann, Hannover

PD Dr. med. Dr. rer. nat. Ingvild Birschmann,
Bad Oeynhausen

Dr. med. Konstantin Kirchmayr, Amberg/Nürnberg

Dr. med. Michael Krause, Leipzig

PD Dr. med. univ. Florian Prüller, Graz

Dr. med. Ute Scholz, Leipzig

We hope this guide helps you introduce integrated, automated, and standardized platelet aggregation testing to your daily laboratory routine.

Abbreviations

Ara	Arachidonic acid
ADP	Adenosine diphosphate
COL	Collagen
COX-1	Cyclooxygenase 1
DNA	Desoxyribonucleic acid
EPI	Epinephrine
GP	Glycoprotein
LTA	Light transmission aggregometry
MPS	Myeloproliferative syndrome
NSAID	Nonsteroidal anti-inflammatory drugs
PAR	Protease-activated receptor
PPP	Platelet-poor plasma
PRP	Platelet-rich plasma
RIS	Ristocetin
RNA	Ribonucleic acid
TRAP	Thrombin-receptor-activating peptide
TxA₂	Thromboxan A ₂
UDM	User-defined method
VWF	von Willebrand factor

Guidelines and literature

Available guidelines

1. Hayward CP, et al. Development of North American consensus guidelines for medical laboratories that perform and interpret platelet function testing using light transmission aggregometry. *Am J Clin Pathol.* 2010;134:955-63. Key aspect: How to interpret LTA results
2. Harrison P, et al. Guidelines for the laboratory investigation of heritable disorders of platelet function. *Br J Haematol.* 2011;155:30-40.
3. Cattaneo M, et al. Recommendations for the standardization of light transmission aggregometry: a consensus of the Working Party from the Platelet Physiology Subcommittee of SSC/ISTH. *J Thromb Haemost.* 2013;11:1183-9. Key aspect: Recommendations for technical procedures
4. Knöfler R, et al. Diagnose von Thrombozytenfunktionsstörungen – Thrombozytopathien GTH e.V.: AWMF-S2K-LL, Register-Nr.: 086–003, Update 2018.
5. Gresele P, for the Subcommittee on Platelet Physiology. Diagnosis of inherited platelet function disorders: guidance from the SSC of the ISTH. *J Thromb Haemost.* 2015;13:314-22.

Additional sources

6. Gawaz M. Das Blutplättchen. Georg Thieme Verlag. ISBN 3-13-118221-0
7. Michelson AD. Platelets. Fourth edition; 2019. ISBN: 978-0-12-813456-6
8. Scharf R. Drugs that affect platelet function. *Semin Thromb Haemost.* 2012;38:865-83.
9. George JM, et al. The clinical importance of acquired abnormalities of platelet function. *NEJM.* 1991;324(1):27-39.

Literature on integrated and automated platelet aggregation

10. Frere C, Kobayashi K, Dunois C, Amiral J, Morange P-E, Alessi M-C. Assessment of platelet function on the routine coagulation analyzer Sysmex CS-2000i. *Platelets.* 2018;29(1):95-7. doi: 10.1080/09537104.2017.1353683 (published online)
11. Bret V-E, Pougault B, Guy A, Castet S, Huguenin Y, Pillois X, James C, Fiore M. Assessment of light transmission aggregometry on the routine coagulation analyzer Sysmex CS-2500 using CE-marked agonists from Hyphen Biomed. *Platelets.* Early online:1-3. doi: 10.1080/09537104.2018.1528346
12. Ling L, Liao J, Niu Q, Wang X, Jia J, Zuo C, Jiang H, Zhou J. Evaluation of an automated light transmission aggregometry. *Platelets.* 28;7:712-9. doi: 10.1080/09537104.2016.1265923 (published online)
13. Platten S, McCormick A, Bukht M, Gurney D, Holding I, Moore GW. A multicenter study to evaluate integrated platelet aggregometry on Sysmex CS-series coagulation analyzers- preliminary findings. *Res Pract Thromb Haemost.* 2018;2:778-89. doi: 10.1002/rth2.12140
14. Stratmann J, Karmal L, Zwinge B, Miesbach W. Platelet aggregation testing on a routine coagulation analyzer: a method comparison study. *Clin Appl Thromb/Hemost.* 2019 Oct;25:1-10. doi: 10.1177/1076029619885184
15. Lawrie AS, Kobayashi K, Lane PJ, Mackie IJ, Machin SJ. The automation of routine light transmission platelet aggregation. *Int J Lab Hemost.* 2014;36:431-8. doi: 10.1111/ijlh.12161
16. Sakayori T, Watanabe Y, Nakajima K, Misawa E, Kobayashi K, Kurono H, Arai N, Takaoka H. Introduction and evaluation of light transmission platelet aggregation method on the Sysmex CS-series Automated Coagulation Analyzer. *Sysmex J Int.* 2016;26(1).

Posters

17. Assessment of PLT Funktion on CS PE Morange France
18. ISTH2015_HYPHEN_PAP8vsCS
SJ Machin_Prec and Throughput CS
19. Prüller F, et al. Platelet function testing using Born's optical aggregometry on automated coagulation analyzer systems compared to a manual aggregometer.

At Siemens Healthineers, we pioneer breakthroughs in healthcare. For everyone. Everywhere. Sustainably. As a leader in medical technology, we want to advance a world in which breakthroughs in healthcare create new possibilities with a minimal impact on our planet. By consistently bringing innovations to the market, we enable healthcare professionals to innovate personalized care, achieve operational excellence, and transform the system of care.

Our portfolio, spanning in vitro and in vivo diagnostics to image-guided therapy and cancer care, is crucial for clinical decision-making and treatment pathways. With the unique combination of our strengths in patient twinning,* precision therapy, as well as digital, data, and artificial intelligence (AI), we are well positioned to take on the greatest challenges in healthcare. We will continue to build on these strengths to help overcome the world's most threatening diseases, enable efficient operations, and expand access to care.

We are a team of more than 71,000 Healthineers in over 70 countries passionately pushing the boundaries of what is possible in healthcare to help improve the lives of people around the world.

**Personalization of diagnosis, therapy selection and monitoring, aftercare, and managing health.*

**This brochure is intended for use outside the U.S. (OUS) only.
Display or promotion to U.S. audiences is prohibited.**

Product availability may vary from country to country and is subject to varying regulatory requirements.

Please contact your local representative for availability.

CN-3000 System refers to Automated Blood Coagulation Analyzer CN-3000. CN-6000 System, CS-2500 System, and CS-5100 System, refer to Automated Blood Coagulation Analyzer CN-6000, CS-2500, and CS-5100, respectively.

Siemens Healthineers Headquarters

Siemens Healthineers AG
Siemensstr. 3
91301 Forchheim, Germany
Phone: +49 9191 18-0
siemens-healthineers.com

Published by

Siemens Healthcare Diagnostics
Products GmbH
Specialty Lab Solutions
Emil-von-Behring-Strasse 76
35041 Marburg, Germany