

**Virtual Touch Tissue Imaging:**

**ACUSON S2000™ ultrasound systems:** There are no prerequisites for first generation Virtual Touch™\* applications. However, release 2.0 software is required for Virtual Touch™ HD applications and to utilize the 9L4 transducer. Other transducers compatible with Virtual Touch are the 4C1 and 4V1.

*\*Not commercially available in the United States*

*Please Note: Additional technical pre-requisites may apply. Upon receiving your request, your local Siemens representative will clarify whether your system meets the requirements.*