syngo Velocity Vector Imaging technology:

ACUSON S2000™ ultrasound system: *syngo*® Velocity Vector Imaging technology requires Cardiac Imaging and Quantification or Cardiac Applications on the ACUSON S2000 system.

ACUSON X300[™] ultrasound system, premium edition (PE): Available on ACUSON X300 ultrasound systems, premium edition only, starting with SN 313001. Also requires any of the Cardiovascular or Advanced Cardiovascular packages available on the ACUSON X300 PE system.

ACUSON and SONOLINE Antares™ systems: Requires Cardiac Applications.

ACUSON Sequoia™ C256 or C512 ultrasound system: Requires Flat Panel Monitor and HSIX on all systems.

ACUSON Sequoia™ 512 system: Requires Cardiac Imaging and Quantification Package as well as Advanced Imaging and Flat Panel Monitor with HSIX.

Please Note: Additional technical pre-requisites may apply. Upon receiving your request, your local Siemens representative will clarify whether your system meets the requirements.