

At Siemens Healthineers, our purpose is to enable healthcare providers to increase value by empowering them on their journey toward expanding precision medicine, transforming care delivery, and improving patient experience, all made possible by digitalizing healthcare.

An estimated 5 million patients globally benefit every day from our innovative technologies and services in the areas of diagnostic and therapeutic imaging, laboratory diagnostics, and molecular medicine, as well as digital health and enterprise services.

We are a leading medical technology company with over 120 years of experience and 18,000 patents globally. Through the dedication of more than 50,000 colleagues in 75 countries, we will continue to innovate and shape the future of healthcare.

Acute Care, ADVIA, ADVIA Centaur, CLINITEK Advantus, CLINITEK Status, DCA Vantage, epoc, RAPIDComm, RAPIDLab, RAPIDPoint, Stratus, Xprecia Stride, and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc., or its affiliates. All other trademarks and brands are the property of their respective owners.

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

Siemens Healthineers Headquarters
Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone: +49 9131 84-0
siemens-healthineers.com

Published by
Siemens Healthcare Diagnostics Inc.
Point of Care Diagnostics
2 Edgewater Drive
Norwood, MA 02062-4637
USA
Phone: +1 781-269-3000

Point of Care Informatics

POCcelerator Device Interface List

siemens-healthineers.com/poc-informatics

Vendor	Device	Connection	Application
4S Information Systems	DAWN ^{NEW}	Unidirectional	Software
77 Elektronika	DocUReader 2 Pro	Bidirectional Download operator/patient Configuration Positive patient identification	Urinalysis
Abaxis	Piccolo Xpress	Unidirectional	Clinical chemistry
Abbott Diabetes Care	FreeStyle Precision Pro (W-LAN)	Bidirectional Download operator/patient Configuration Positive patient identification	Glucose Ketone
	PCx Precision	Unidirectional	Glucose
	Precision Xceed Pro	Bidirectional Download operator/patient Configuration Positive patient identification	Glucose Ketone
Abbott Point of Care	Afinion AS 100 Analyzer	Bidirectional Configuration Download operator	HbA1c
	Afinion 2	Bidirectional Configuration Download operator	HbA1c
	Alere i	Unidirectional	Multitest
	ID NOW	Bidirectional Download operator Configuration Positive Patient identification	Multitest
	CELL-DYN Emerald	Unidirectional	Hematology
	Cholestech LDX System	Unidirectional	Multitest
	DIGIVAL (Alere Reader)	Unidirectional	Multitest
	Pima	Unidirectional	HIV
	i-STAT 1	Bidirectional (via i-STAT-DE-DMI software) Download operator/patient Configuration Positive patient identification	Multitest
	i-STAT Alinity	Bidirectional Download operator/patient Configuration Positive patient identification	Multitest
i-STAT/DE/DMI	Bidirectional Download operator/patient Configuration Positive patient identification	Supplier software	

Vendor	Device	Connection	Application
Acon Laboratories	Mission U120	Unidirectional	Urinalysis
Aidian (former Orion Diagnostica)	QuikRead 101	Unidirectional	CRP
	QuikRead go	Bidirectional Download operator/patient Configuration	CRP
A.MENARINI Diagnostics	Zenit RA	Unidirectional	Immunoassay
Analyticon	CombiScan 100	Unidirectional	Urinalysis
	Urilyzer 100	Bidirectional Download operator	Urinalysis
	Urilyzer 100 Pro	Bidirectional Configuration Download operator/patient Positive patient identification	Urinalysis
Arkray	ADAMS HA-8160	Unidirectional	HbA1c
Beckman Coulter	Coulter Ac-T diff	Unidirectional	Hematology (no QCI)
	Coulter Ac-T diff2	Unidirectional	Hematology
	Coulter Ac-T 8diff	Unidirectional	Hematology
	Coulter Ac-T 10diff	Unidirectional	Hematology
	DxH 500	Unidirectional	Hematology
Becton Dickinson	BD Veritor Plus	Unidirectional (via BD Cloud Connectivity Solution)	Multitest
biosurfit	spinit	Unidirectional	Multitest
Boditech	AFIAS-1	Unidirectional	Immunoassay
	AFIAS-6	Unidirectional	Immunoassay
Boule Medical	Medonic M-series M16/M20	Unidirectional	Hematology
	Medonic M-series M16M/M20M	Unidirectional	Hematology
	Medonic M-series M16C/M20C	Unidirectional	Hematology
	Medonic M-series M16S/M20S	Unidirectional	Hematology
	Medonic M-series M32B	Unidirectional	Hematology
	Medonic M-series M32M	Unidirectional	Hematology
	Medonic M-series M32C	Unidirectional	Hematology
	Medonic M-series M32S	Unidirectional	Hematology
	Swelab Alpha / Alpha Plus Basic	Unidirectional	Hematology
	Swelab Alpha / Alpha Plus Standard	Unidirectional	Hematology
	Swelab Alpha / Alpha Plus Cap	Unidirectional	Hematology
Swelab Alpha / Alpha Plus Sampler	Unidirectional	Hematology	
Caretaker Medical	Caretaker4 ^{NEW}	Unidirectional	Vital Signs Monitor
Cepheid	GeneXpert ^{NEW}	Unidirectional	Multitest
	GeneXpert Infinity ^{NEW}	Unidirectional	Multitest

Vendor	Device	Connection	Application
Cepheid	GeneXpert Xpress ^{NEW}	Bidirectional Configuration Download operator	Multitest
Chirus	DXpress Reader	Unidirectional	Multitest
concile	Q100	Unidirectional	Immunoassay
DnaNudge	NudgeBox ^{NEW}	Unidirectional	Immunoassay
Dräger	Jaundice Meter JM-105	Unidirectional	Bilirubin
Drucker Diagnostics	QBC STAR	Unidirectional	Hematology
EKF Diagnostics	Quo-Lab via CIB ^{NEW}	Bidirectional Download operator/patient Configuration	HbA1c
	Quo-Test ^{NEW}	Bidirectional	HbA1c
Enigma Diagnostics	Enigma MiniLab	Unidirectional	Multitest
ET Healthcare	ET Pylon ^{NEW}	Unidirectional	Multitest
Eurolyser Diagnostica	Eurolyser cube	Unidirectional	Multitest
	Eurolyser cube S	Unidirectional	Multitest
	Eurolyser Solo	Unidirectional	Multitest
	Eurolyser Smart	Unidirectional	Multitest
Fuji	FUJI DRI-CHEM 4000i	Unidirectional	Clinical chemistry
	FUJI DRI-CHEM NX500	Unidirectional	Clinical chemistry
	FUJI DRI-CHEM NX10N	Unidirectional	Clinical chemistry
	FUJI DRI-CHEM NX700	Unidirectional	Clinical chemistry
GenePOC Diagnostics	Revogene	Bidirectional Download operator Configuration	Multitest
Haemonetics	TEG 6s	Bidirectional Positive patient identification	Thromlastometry
Helena Laboratories	Cascade Abrazo	Bidirectional Configuration	Coagulation
	Cascade POC	Bidirectional Download operator	Coagulation
HemoCue	Glucose 201 DM	Bidirectional Download operator/patient Configuration Positive patient identification	Glucose
	Glucose 201 DM RT	Bidirectional Download operator/patient Configuration Positive patient identification	Glucose
	HbA1c 501	Unidirectional	Hemoglobin
	Hb 201 DM	Bidirectional Download operator/patient Configuration Positive patient identification	Hemoglobin
	WBC Diff	Unidirectional	Hematology
HemoSonics	Quantra	Bidirectional Download operator/patient	Coagulation
HLS Therapeutics Inc.	CSAN Pronto ^{NEW}	Unidirectional	Hematology
Hologic	Rapid fFN	Unidirectional	Immunoassay

Vendor	Device	Connection	Application
Hologic	PeriLynx	Unidirectional	Immunoassay
HORIBA	ABX Micros 60	Unidirectional	Hematology
	ABX Micros ES 60	Unidirectional	Hematology
	ABX Microsemi CRP	Unidirectional	Multitest
	Yumizen H500	Unidirectional	Hematology
Lite-On Technology Corp.	skyla HB1	Bidirectional Download operator	Clinical chemistry
	skyla Hi	Bidirectional Download operator	HbA1c
LumiraDx	LumiraDx ^{NEW}	Bidirectional Download operator Positive patient identification Configuration	Multitest
Macherey-Nagel	URYXXON Relax (Connect)	Bidirectional Positive patient identification	Urinalysis
Medtronic	ACT Plus	Unidirectional	Coagulation
	HMS Plus	Unidirectional	Coagulation
Mitsubishi	PATHFAST	Unidirectional	Immunoassay
Nal von minden	Rapid Slide Scanner III ^{NEW} (RSS)	Unidirectional	Drug-Screen
Nihon Kodan	Celltac Alpha (MEK 6500/6510)	Unidirectional	Hematology
	Celltac Chemi (CHM-4100K)	Unidirectional	Hematology
Norma Diagnostika	Norma Icon 3	Unidirectional	Hematology
	Norma Icon 5	Unidirectional	Hematology
Nova Biomedical	NovaNet*	Bidirectional Download operator Positive patient identification	Supplier software
	Stat Profile pHox/pHOx Plus	Bidirectional Positive patient identification	Blood gas
	Stat Profile PHox Ultra	Unidirectional	Blood gas
	Stat Profile CCX	Bidirectional Positive patient identification	Blood gas
	Stat Profile Prime	Unidirectional	Blood gas
	Stat Profile Prime Plus	Bidirectional Download operator Positive patient identification (via NovaNet)	Blood gas
	StatStrip Glucose	Bidirectional Download operator/patient Configuration Positive patient identification	Glucose
	StatStrip Glucose/Ketone	Bidirectional Download operator/patient Configuration Positive patient identification	Glucose Ketone

Vendor	Device	Connection	Application
Nova Biomedical	StatStrip Hb/Hct	Bidirectional Download operator Positive patient identification (via NovaNet)	Hemoglobin
	StatStrip HHL ^{NEW}	Bidirectional Download operator Positive patient identification (via NovaNet)	Hemoglobin/ Hematocrit/Lactate
	StatStrip Lactate	Bidirectional Download operator/patient Configuration	Lactate
	StatSensor Creatinine	Bidirectional Download operator/patient Configuration	Creatinine
OPTI Medical	OPTI CCA-TS2	Bidirectional Download operator/patient	Blood gas
Philips	BiliChek	Unidirectional	Bilirubin
PixCell Medical Technologies	Hemoscreen	Unidirectional	Hematology
Quidel	Sofia	Unidirectional	Immunoassay
	Sofia 2 ^{NEW}	Unidirectional	Immunoassay
	Triage MeterPlus	Bidirectional Download operator	Cardiac markers
	Triage MeterPro	Bidirectional Download operator	Cardiac markers
Radiometer	ABL7xx series	Bidirectional Positive patient identification	Blood gas
	ABL8xx series	Bidirectional Positive patient identification	Blood gas
	ABL80 FLEX series	Bidirectional Positive patient identification	Blood gas
	ABL90 FLEX/FLEX Plus	Bidirectional Positive patient identification	Blood gas
	AQT90 FLEX	Bidirectional Positive patient identification	Multitest
	AQURE	Bidirectional Positive patient identification Download operator	Supplier software
	NPT7	Bidirectional Positive patient identification	Blood gas
	RADIANCE	Bidirectional Positive patient identification Download operator	Supplier software
	Response Biomedical	RAMP 200	Bidirectional Download operator Configuration
Roche Diagnostics	ACCU-CHEK Inform II	Bidirectional Download operator/patient Configuration Positive patient identification	Glucose

*NovaNet is needed for the StatStrip-Family connections with a firmware version greater than 4.17.8.47.

Vendor	Device	Connection	Application	
Roche Diagnostics	CARDIAC reader	Bidirectional Configuration	Cardiac markers	
	CoaguChek Pro II	Bidirectional Download operator/patient Configuration Positive patient identification	Coagulation	
	CoaguChek XS	Unidirectional	Coagulation	
	CoaguChek XS Plus	Bidirectional Download operator/patient Configuration Positive patient identification	Coagulation	
	CoaguChek XS Pro	Bidirectional Download operator/patient Configuration Positive patient identification	Coagulation	
	cobas b101	Bidirectional Download operator/patient Configuration Positive patient identification	Hemoglobin	
	cobas b 123	Bidirectional Download operator/patient Positive patient identification	Blood gas	
	cobas c 111	Unidirectional	HbA1c	
	cobas h 232	Bidirectional Download operator/patient Positive patient identification Configuration	Cardiac markers	
	cobas Liat	Bidirectional Download operator	Multitest	
	cobas 6000	Unidirectional	Clinical chemistry	
	cobas IT 1000	Bidirectional Download operator/patient Positive patient identification	Supplier software	
	Reflotron Plus	Unidirectional	Clinical chemistry	
	Reflotron Sprint	Unidirectional	Clinical chemistry	
	Urisys 1100	Bidirectional (with firmware version 5.x or higher) Download operator	Urinalysis	
	Samsung	LABGEO ^{HC10}	Bidirectional Download operator	Hematology
		LABGEO ^{IB10}	Bidirectional Download operator	Immunoassay
LABGEO ^{PT10}		Bidirectional Download operator	Clinical chemistry	
SD Biosensor	Standard F200	Unidirectional	Multitest	
	Standard F2400	Unidirectional	Multitest	
Siemens Healthcare Diagnostics	ADVIA [®] 360 Hematology System	Unidirectional	Hematology	
	ADVIA Centaur [®] CP/XP Immunoassay System	Unidirectional	Immunoassay	
	CLINITEK Advantus [®] Urine Chemistry Analyzer	Unidirectional	Urinalysis	

Vendor	Device	Connection	Application	
Siemens Healthcare Diagnostics	CLINITEK Status [®] Analyzer	Unidirectional	Urinalysis	
	CLINITEK Status [®] + Urine Chemistry Analyzer	Bidirectional Download operator Configuration	Urinalysis	
	DCA Vantage [®] Analyzer	Bidirectional Download operator	HbA1c	
	epoc [®] Blood Analysis System	Bidirectional Download operator Configuration Positive patient identification	Blood gas	
	RAPIDComm [®] Data Management System	Bidirectional Download operator Positive patient identification	Supplier software	
	Atellica 1500 [®] Automated Urinalysis System	Bidirectional Configuration Positive sample identification	Urinalysis	
	Atellica [®] UAS 800 Automated Urinalysis System	Bidirectional Configuration Positive sample identification	Urinalysis	
	RAPIDLab [®] 348EX Blood Gas Systems	Unidirectional	Blood gas	
	RAPIDLab 1200 series Blood Gas Systems	Bidirectional Download operator Configuration Positive patient identification	Blood gas	
	RapidPoint 400, 405 Blood Gas Systems	Bidirectional Download operator Configuration Positive patient identification	Blood gas	
	RAPIDPoint 500 Blood Gas Systems	Bidirectional Download operator Configuration Positive patient identification	Blood gas	
	RAPIDPoint 500e Blood Gas Systems	Bidirectional Download operator Configuration Positive patient identification	Blood gas	
	Stratus [®] CS Acute Care [™] Diagnostic System	Unidirectional	Cardiac markers	
	Stratus CS 200 Acute Care Diagnostic System	Bidirectional Download operator	Cardiac markers	
	Xprecia Stride [™] Coagulation Analyzer	Bidirectional Download operator Configuration	Coagulation	
	Sight Diagnostics	OLO	Unidirectional	Hematology
	Sphere Medical	Proxima	Unidirectional	Blood gas
Stat-Diagnostica	DiagCORE	Unidirectional	Multitest	
Sysmex	Alifax Roller	Unidirectional	Blood gas	
	CA-560	Unidirectional	Coagulation	
	CA-600	Unidirectional	Coagulation	
	pocH-100i	Unidirectional	Hematology	

Vendor	Device	Connection	Application
Sysmex	XN-350, XN-450, XN-550	Unidirectional (via SUIT Manager)	Hematology
	XP-300	Unidirectional	Hematology
	XS-800i	Unidirectional (via SUIT Manager)	Hematology
Thermo Fisher Scientific	B·R·A·H·M·S Direct Reader	Bidirectional Configuration	Multitest
Veda.Lab	Easy Reader Plus	Bidirectional Download operator	Multitest
Welch Allyn	Connex Spot Monitor ^{NEW}	Unidirectional	Vital Signs Monitor
Werfen	AVOXimeter 1000e	Unidirectional	Blood gas
Werfen	AVOXimeter 4000	Unidirectional	Blood gas
	BIO-FLASH	Unidirectional	Immunoassay
	GEM Premier 3000	Bidirectional Positive patient identification	Blood gas
	GEM Premier 3500	Bidirectional Positive patient identification	Blood gas
	GEM Premier 4000/5000	Bidirectional Download operator Positive patient identification	Blood gas
	GEMweb Plus	Bidirectional Download operator Positive patient identification	Supplier software
	HEMOCHRON Jr. Signature	Unidirectional	Coagulation
	HEMOCHRON Signature Elite	Bidirectional Download operator Configuration	Coagulation
	HEMOCHRON Response	Unidirectional	Coagulation
	ProTime InRhythm	Unidirectional	Coagulation
	ROTEM delta	Bidirectional Positive patient identification	Thromlastometry
	ROTEM gamma	Bidirectional Positive patient identification	Thromlastometry
	ROTEM sigma	Bidirectional Positive patient identification	Thromlastometry